






Primary analysis metrics from the 1 x 150 bp and 2 x 75 bp cycle sequencing runs (Table 1) demonstrate the trade-offs between the 2 run configurations. The overall quality of the 1 x 150 bp run is higher with 94% of bases greater than or equal to Q30, compared to 86% of bases greater than or equal to Q30 for the 2 x 75 bp runs. This is also reflected in the average PhiX error rate of 0.33% for 1 x 150 bp compared to 0.45% average with 2 x 75 bp runs. While performance of the single-read runs is higher than the paired-end runs for these metrics, paired-end runs have the added benefit of enabling identification of duplicates that can be removed in subsequent analyses.



**Figure 4:** Box plot over all tiles of percentage of base calls that are errors in spiked-in PhiX in each cycle—All clusters passing filters are included from reads assigned to spiked-in control DNA. (A) 1 x 150 bp single-end run and (B) 2 x 75 bp paired-end runs performed on the HiSeq 4000 System.

## Secondary Data Analysis

Multiple secondary analyses were performed using the MethylSeq v1.0 BaseSpace® App.<sup>11</sup> From each run, data from multiple lanes were pooled and analyzed to demonstrate the number of samples that can be processed per flow cell according to the desired depth of coverage.

Overall, the coverage obtained using the 1 x 150 bp sequencing run is higher than the 2 x 75 bp run for an equivalent number of lanes (Table 2). There are 2 main reasons for the increased coverage. The first is due to duplicates, which are not removed in the single-read 1 x 150 bp sequencing run. In the paired-end run, duplicate reads are easily distinguished and can be removed from subsequent analyses. The second reason for the lower coverage of the 2 x 75 bp run is the slight reduction in data quality at the end of Read 2, which yields a lower percentage of aligned reads. The other metrics for percent methylation show a high degree of reproducibility between analyses from different run configurations and from combining different numbers of lanes.

**Table 2: Primary Analysis Sequencing Metrics from 1 x 150 bp and 2 x 75 bp cycle WGBS Sequencing Runs**

Read Length (bp)	1 x 150 (2 lanes)	2 x 75 (2 lanes)	1 x 150 (3 lanes)	2 x 75 (3 lanes)
% Duplicates	N/A	21.0	N/A	24.4
% Aligned Reads	63.8	59.2	64.6	58.2
Coverage	21.7	15.7	32.5	22.7
% Methylated Cs in CpG	51.4	51.5	51.4	51.5
% Methylated Cs in CHG	0.7	0.4	0.7	0.4
% Methylated Cs in CHH	0.5	0.5	0.5	0.5

## Guidelines

WGBS involves 2 significant challenges that result from the use of bisulfite converted DNA. First, the conversion of cytosines to thymines causes a reduction in genome diversity that impacts the quality of sequencing data. Second, the conditions required for complete bisulfite conversion result in significant DNA degradation. The first challenge is addressed with the “spike-in” addition of DNA from a well balanced genome, enabling accurate base calling on the sequencer. The TruSeq DNA Methylation Kit addresses the second challenge by capturing ssDNA fragments and reducing sample loss. Using a high-quality library from this kit, while following Illumina-recommended parameters, users can obtain highly accurate MethylSeq results with the HiSeq 3000/HiSeq 4000 Systems.

The HiSeq 3000/HiSeq 4000 Systems provide the option for single-read or paired-end runs. Due to the average size of library products from bisulfite converted DNA, read configurations are optimally set at either 1 x 150 bp or 2 x 75 bp. Both configurations were compared for accuracy as indicated by Q-scores (Figure 3 and Table 1). Low Q-scores decrease the percentage of bases aligned, and increase the amount of data being filtered in the secondary analysis pipeline. Although both run configurations provide high-quality data, the 1 x 150 bp configuration yields slightly higher data accuracy, and increased coverage. However, paired-end runs enable the removal of duplicate reads that can affect subsequent analyses.

