

TruSeq Nano DNA Library Prep for NeoPrep

Protocol Guide

For Research Use Only. Not for use in diagnostic procedures.

Fragment DNA	3
Prepare Samples for Loading	4
Set Up Run and Load Library Card	5
Unload Libraries	7
[Optional] Validate Libraries Manually	8
[Optional] Normalize Libraries Manually	9
Pool Libraries	10
Acronyms	11
Technical Assistance	13

ILLUMINA PROPRIETARY

Document # 15059579 v01
October 2015

illumina®



This document and its contents are proprietary to Illumina, Inc. and its affiliates ("Illumina"), and are intended solely for the contractual use of its customer in connection with the use of the product(s) described herein and for no other purpose. This document and its contents shall not be used or distributed for any other purpose and/or otherwise communicated, disclosed, or reproduced in any way whatsoever without the prior written consent of Illumina. Illumina does not convey any license under its patent, trademark, copyright, or common-law rights nor similar rights of any third parties by this document.

The instructions in this document must be strictly and explicitly followed by qualified and properly trained personnel in order to ensure the proper and safe use of the product(s) described herein. All of the contents of this document must be fully read and understood prior to using such product(s).

FAILURE TO COMPLETELY READ AND EXPLICITLY FOLLOW ALL OF THE INSTRUCTIONS CONTAINED HEREIN MAY RESULT IN DAMAGE TO THE PRODUCT(S), INJURY TO PERSONS, INCLUDING TO USERS OR OTHERS, AND DAMAGE TO OTHER PROPERTY.

ILLUMINA DOES NOT ASSUME ANY LIABILITY ARISING OUT OF THE IMPROPER USE OF THE PRODUCT(S) DESCRIBED HEREIN (INCLUDING PARTS THEREOF OR SOFTWARE).

© 2015 Illumina, Inc. All rights reserved.

Illumina, 24sure, BaseSpace, BeadArray, BlueFish, BlueFuse, BlueGnome, cBot, CSPPro, CytoChip, DesignStudio, Epicentre, ForenSeq, Genetic Energy, GenomeStudio, GoldenGate, HiScan, HiSeq, HiSeq X, Infinium, iScan, iSelect, MiSeq, MiSeqDx, MiSeq FGx, NeoPrep, NextBio, Nextera, NextSeq, Powered by Illumina, SureMDA, TruGenome, TruSeq, TruSight, Understand Your Genome, UYG, VeraCode, verifi, VeriSeq, the pumpkin orange color, and the streaming bases design are trademarks of Illumina, Inc. and/or its affiliate(s) in the U.S. and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.

AMPure®, Beckman, and Beckman Coulter are trademarks or registered trademarks of Beckman Coulter, Inc.

Fragment DNA

Preparation

- 1 Turn on the Covaris instrument and follow manufacturer instructions to set up your instrument.

Procedure

- 1 Quantify gDNA samples using a fluorometric-based method that uses dsDNA binding dyes such as Qubit or QuantiFlour.
- 2 Normalize gDNA with RSB in a final volume of 15 μ l in separate wells of a new PCR plate:
 - ▶ For a 350 bp insert size—1.67 ng/ μ l of sample
 - ▶ For a 550 bp insert size—5 ng/ μ l of sample
- 3 Transfer 15 μ l of each normalized DNA to a microTUBE-15.
- 4 Centrifuge, using a microTUBE adapter, at $3000 \times g$ for 1 minute.
- 5 Fragment the DNA on a Covaris.

Table 1 Covaris S220 or E220 Settings

Setting	350 bp Insert	550 bp Insert
Duty factor	20%	
Peak Incident Power	18 W	
Cycles per burst	50	
Duration	45 seconds	22 seconds
Temperature	20°C	
Water Level—S220	15	
Water Level—E220	10	

Table 2 Covaris M220 Settings

Setting	350 bp Insert	550 bp Insert
Duty factor	20%	
Peak Incident Power	30 W	
Cycles per burst	50	
Duration	42 seconds	23 seconds
Temperature	20°C	

- 6 Centrifuge, using a microTUBE adapter, at $600 \times g$ for 5 seconds.
- 7 Transfer 15 μ l fragmented DNA from each microTUBE-15 to a separate well of a new PCR plate.

SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to 7 days.

Prepare Samples for Loading

Procedure

- 1 Add 1 of the following to a new 1.5 ml microcentrifuge tube:
 - ▶ For a 350 bp insert size—900 μ l SC350
 - ▶ For a 550 bp insert size—700 μ l SC550
- 2 Vortex DMB until well-dispersed.
- 3 Add 100 μ l DMB to the microcentrifuge tube containing SC350 or SC550. Vortex for 5 seconds.
- 4 Pour the DMB and SC mixture into a new reagent reservoir.
- 5 Add 35 μ l DMB and SC mixture to each well of the sample plate. Pipette to mix.
- 6 Shake or vortex at 1400 rpm for 12 minutes.

Set Up Run and Load Library Card

Procedure

- 1 Vortex the reagent plate for 3 seconds.
- 2 Centrifuge at $600 \times g$ for 5 seconds.
- 3 Select **Prepare Libraries** on the NeoPrep System Welcome screen.
- 4 Do the following and then select **Next**.
 - ▶ If running in BaseSpace mode, select a run.
 - ▶ If running in standalone mode, use the following options to select a protocol:
 - ▶ Select **Select by barcode**, and then scan the reagent plate barcode or enter the reagent plate serial number.
 - ▶ Select **Select by name**, and then select **TruSeq Nano DNA**.
- 5 Configure the run. Select **Next**.
 - ▶ Select the Insert Size.
- 6 Review the run and sample information. Select **Next**.
- 7 Enter the consumable tracking information. Select **Next**.
- 8 Place the library card on the library card stage.



WARNING

To avoid instrument damage, make sure that the library card guide is not on the library card.

- 9 Close the library card compartment door. Select **Verify Library Card**.
- 10 Place the library card guide on the library card.
- 11 Load the entire contents of the oil vial into the library card using the oil funnel.



WARNING

Use the required pipette tips. Other tips are not supported and can result in reagents not dispensing properly and run failure.

The loading angle of the pipette depends on the item being dispensed. The angle is specified in each step of the NPCS loading guide and is depicted in these procedures.

- 12 Transfer 45 μ l of prepared samples 1–8.
- 13 Transfer 45 μ l of prepared samples 9–16.
- 14 If you are preparing < 16 samples, add 45 μ l RSB to empty sample wells.
- 15 Transfer 125 μ l of the large reagents i–iv.
- 16 Transfer 125 μ l of the large reagents v–vii.
- 17 Vortex DMB until well-dispersed.
- 18 Add 60 μ l DMB to the large reagent well viii.
- 19 Transfer 15 μ l of small reagents 1–4, and then 5–8.
- 20 Transfer 5 μ l of small reagents a–d, and then e–h.
- 21 Transfer 3 μ l of adapters A–H.
- 22 Transfer 3 μ l of adapters I–P.

23 Remove the library card guide.



WARNING

To avoid instrument damage, make sure that the library card guide is removed from the library card.

24 Close the library card compartment door. Select **Start Run**.

25 When the run is complete, select **Next**.

Unload Libraries



WARNING

The used library card contains hazardous materials. Personal injury can occur through inhalation, ingestion, skin contact, and eye contact. Wear protective equipment, including eye protection, gloves, and a laboratory coat. Handle the used library card as chemical waste. Dispose of containers and any unused contents in accordance with the governmental safety standards for your region. For more information, see the SDS for this kit at support.illumina.com/sds.html.

Procedure

- 1 Add 10 μ L RSB to each well of a new PCR plate labeled 1–16.
- 2 Open the library card compartment door and place the library card guide on the library card.
- 3 Use a 200 μ L pipette to transfer 20 μ L from library card collection wells 1L–8L, and then 9L–16L to corresponding wells 1–16 of the plate. Pipette to mix.
- 4 Centrifuge briefly.
- 5 Transfer the entire volume from plate wells 1–8, and then 9–16 to the center indent in the membrane of the corresponding library separation tubes 1–16.
- 6 Let stand for 10 seconds while the oil is absorbed in the tubes.
- 7 Transfer the entire volume from library separation tubes 1–8, and then 9–16 to the corresponding wells 1–16 of a new PCR plate.
- 8 Remove the library card and library card guide from the library card stage.
- 9 Discard the library card in accordance with applicable standards.
- 10 Close the library card compartment door, and then select **Home**.
- 11 Select from the following options:

Table 3 Post Run Options

NeoPrep System Quantification	NeoPrep System Normalization	Pooling Required	Then...
Yes	Yes	No	The protocol stops here. The final library is normalized to 10 nM. Proceed to cluster generation.
Yes	Yes	Yes	Proceed to <i>Pool Libraries</i> .
Yes	No	Yes or No	Proceed to <i>[Optional] Normalize Libraries Manually</i> .
No	No	Yes or No	Proceed to <i>[Optional] Validate Libraries Manually</i> .

SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to 2 months.

[Optional] Validate Libraries Manually

Procedure

- 1 Quantify the libraries using a fluorometric quantification method that uses dsDNA binding dyes or qPCR.
- 2 Check the library size distribution:
 - ▶ If using a High Sensitivity DNA chip:
 - ▶ Dilute the DNA library 1:10 with water.
 - ▶ Run 1 μ l diluted DNA library.
 - ▶ If using a DNA 7500 chip, run 1 μ l undiluted DNA library.

SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to 2 months.

[Optional] Normalize Libraries Manually

Procedure

- 1 Transfer 5 μ l from each well of the library plate to the corresponding wells of a midi plate.
- 2 Normalize each library to 10 nM with Tris-HCl 10 mM, pH 8.5 with 0.1% Tween 20. Pipette to mix.
- 3 Select from the following options:
 - ▶ For libraries that do not require pooling, the protocol stops here. Proceed to cluster generation.
 - ▶ For libraries that require pooling, proceed to *Pool Libraries*.

SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to 2 months.

Pool Libraries

Procedure

- 1 Determine the number of samples to combine for each pool.
- 2 Transfer 5 μ l of each library to be pooled from the library plate to a single well of a new PCR plate. Pipette to mix.
- 3 Proceed to cluster generation.

SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to 2 months.

Acronyms

Acronym	Definition
DMB	Digital Microfluidics Beads
RSB	Resuspension Buffer
SC350	Sample Concentration Solution 350
SC550	Sample Concentration Solution 550

Notes

Technical Assistance

For technical assistance, contact Illumina Technical Support.

Table 4 Illumina General Contact Information

Website	www.illumina.com
Email	techsupport@illumina.com

Table 5 Illumina Customer Support Telephone Numbers

Region	Contact Number	Region	Contact Number
North America	1.800.809.4566	Italy	800.874909
Australia	1.800.775.688	Netherlands	0800.0223859
Austria	0800.296575	New Zealand	0800.451.650
Belgium	0800.81102	Norway	800.16836
Denmark	80882346	Spain	900.812168
Finland	0800.918363	Sweden	020790181
France	0800.911850	Switzerland	0800.563118
Germany	0800.180.8994	United Kingdom	0800.917.0041
Ireland	1.800.812949	Other countries	+44.1799.534000

Safety data sheets (SDSs)—Available on the Illumina website at support.illumina.com/sds.html.

Product documentation—Available for download in PDF from the Illumina website. Go to support.illumina.com, select a product, then select **Documentation & Literature**.



Illumina

San Diego, California 92122 U.S.A.

+1.800.809.ILMN (4566)

+1.858.202.4566 (outside North America)

techsupport@illumina.com

www.illumina.com