


TruSeq Small RNA Library Prep

Protocol Guide

For Research Use Only. Not for use in diagnostic procedures.

| | |
|--|----|
| Ligate Adapters | 3 |
| Reverse Transcribe and Amplify Libraries | 5 |
| Purify cDNA Construct | 7 |
| Check Libraries | 9 |
| Normalize Libraries | 10 |
| Acronyms | 11 |
| Technical Assistance | |



ILLUMINA PROPRIETARY

illumina®

This document and its contents are proprietary to Illumina, Inc. and its affiliates ("Illumina"), and are intended solely for the contractual use of its customer in connection with the use of the product(s) described herein and for no other purpose. This document and its contents shall not be used or distributed for any other purpose and/or otherwise communicated, disclosed, or reproduced in any way whatsoever without the prior written consent of Illumina. Illumina does not convey any license under its patent, trademark, copyright, or common-law rights nor similar rights of any third parties by this document.

The instructions in this document must be strictly and explicitly followed by qualified and properly trained personnel in order to ensure the proper and safe use of the product(s) described herein. All of the contents of this document must be fully read and understood prior to using such product(s).

FAILURE TO COMPLETELY READ AND EXPLICITLY FOLLOW ALL OF THE INSTRUCTIONS CONTAINED HEREIN MAY RESULT IN DAMAGE TO THE PRODUCT(S), INJURY TO PERSONS, INCLUDING TO USERS OR OTHERS, AND DAMAGE TO OTHER PROPERTY.

ILLUMINA DOES NOT ASSUME ANY LIABILITY ARISING OUT OF THE IMPROPER USE OF THE PRODUCT(S) DESCRIBED HEREIN (INCLUDING PARTS THEREOF OR SOFTWARE).

© 2016 Illumina, Inc. All rights reserved.

Illumina, TruSeq, the pumpkin orange color, and the streaming bases design are trademarks of Illumina, Inc. and/or its affiliate(s) in the U.S. and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.

Ligate Adapters

Preparation

- 1 Preheat a thermal cycler to 70°C.
- 2 Choose the thermal cycler preheat lid option and set to 100°C.

Procedure

- 1 Combine the following volumes in a new 200 µl PCR tube on ice:
 - ▶ RA3 (1 µl)
 - ▶ 1 µg total RNA in nuclease-free water (5 µl)
- 2 Pipette to mix, and then centrifuge briefly.
- 3 Place on the preheated thermal cycler.
- 4 Incubate at 70°C for 2 minutes.
- 5 Remove from the thermal cycler and place on ice.
- 6 Preheat the thermal cycler to 28°C.
- 7 Combine the following volumes in a new 200 µl PCR tube on ice. Multiply each volume by the number of samples being prepared. Make 10% extra reagent if you are preparing multiple samples.
 - ▶ HML (2 µl)
 - ▶ RNase Inhibitor (1 µl)
 - ▶ T4 RNA Ligase 2, Deletion Mutant (1 µl)
- 8 Pipette to mix, and then centrifuge briefly.
- 9 Add 4 µl to the tube of RA3/total RNA mixture.
- 10 Pipette to mix.
- 11 Place on the preheated thermal cycler.
- 12 Incubate at 28°C for 1 hour.
- 13 Add 1 µl STP and pipette to mix.
- 14 Continue incubating at 28°C for 15 minutes.
- 15 Remove from the thermal cycler and place on ice.
- 16 Preheat the thermal cycler to 70°C.
- 17 Add $1.1 \times N$ µl RA5 to a new 200 µl PCR tube.
- 18 Place on the preheated thermal cycler.
- 19 Incubate at 70°C for 2 minutes.
- 20 Remove from the thermal cycler and place on ice.
- 21 Preheat the thermal cycler to 28°C.
- 22 Add $1.1 \times N$ µl 10mM ATP to the tube of RA5.
- 23 Pipette to mix.
- 24 Add $1.1 \times N$ µl T4 RNA Ligase to the RA5/ATP mixture.

- 25 Pipette to mix.
- 26 Add 3 μ l to the tube of RA3 mixture.
- 27 Pipette to mix.
- 28 Place on the preheated thermal cycler.
- 29 Incubate at 28°C for 1 hour.
- 30 Remove from the thermal cycler and place on ice.

Reverse Transcribe and Amplify Libraries

Preparation

- 1 Preheat the thermal cycler to 70°C.
- 2 Choose the thermal cycler preheat lid option and set to 100°C.
- 3 Label a new 200 µl PCR tube 12.5 mM dNTP Mix.

Procedure

- 1 Combine the following volumes in the 12.5 mM dNTP Mix tube to dilute to 12.5 mM. Multiply each volume by the number of samples being prepared. Prepare 10% extra reagent if you are preparing multiple libraries.
 - ▶ 25 mM dNTP Mix (0.5 µl)
 - ▶ Ultrapure water (0.5 µl)
- 2 Pipette to mix, and then centrifuge briefly.
- 3 Set aside on ice.
- 4 Add 6 µl each adapter-ligated RNA library to a new 200 µl PCR tube.
- 5 Add 1 µl RNA RT Primer to the tube of adapter-ligated RNA.
- 6 Pipette to mix, and then centrifuge briefly.
- 7 Place on the preheated thermal cycler.
- 8 Incubate at 70°C for 2 minutes.
- 9 Remove from the thermal cycler and place on ice.
- 10 Preheat the thermal cycler to 50°C.
- 11 Combine the following volumes in a new 200 µl PCR tube on ice. Multiply each volume by the number of libraries being prepared. Make 10% extra reagent if you are preparing multiple libraries.
 - ▶ 5X First Strand Buffer (2 µl)
 - ▶ 12.5 mM dNTP Mix (0.5 µl)
 - ▶ 100 mM DTT (1 µl)
 - ▶ RNase Inhibitor (1 µl)
 - ▶ SuperScript II Reverse Transcriptase (1 µl)
- 12 Pipette to mix, and then centrifuge briefly.
- 13 Add 5.5 µl to the tube of adapter-ligated RNA/primer mix.
- 14 Pipette to mix, and then centrifuge briefly.
- 15 Place on the preheated thermal cycler.
- 16 Incubate at 50°C for 1 hour.
- 17 Remove from the thermal cycler and place on ice.
- 18 Combine the following reagents in a new 200 µl PCR tube on ice to prepare the PCR master mix. Multiply each volume by the number of libraries being prepared. Make 10% extra reagent if you are preparing multiple libraries with the same index.
 - ▶ Ultrapure water (8.5 µl)
 - ▶ PML (25 µl)

- ▶ RP1 (2 µl)
 - ▶ RPIX (2 µl)
- 19 Pipette to mix, and then centrifuge briefly.
 - 20 Place on ice.
 - 21 Add 37.5 µl PCR master mix to the adapter-ligated RNA mixture.
 - 22 Pipette to mix, and then centrifuge briefly.
 - 23 Place on ice.
 - 24 Place on the preheated thermal cycler.
 - 25 Incubate using the following program on the thermal cycler:
 - ▶ Choose the preheat lid option and set to 100°C.
 - ▶ 98°C for 30 seconds
 - ▶ 11 cycles of:
 - ▶ 98°C for 10 seconds
 - ▶ 60°C for 30 seconds
 - ▶ 72°C for 15 seconds
 - ▶ 72°C for 10 minutes
 - ▶ 4°C hold
 - 26 Run each library on a High Sensitivity DNA chip.

SAFE STOPPING POINT

If you are stopping, cap the tube and store at -25°C to -15°C for up to 7 days.

Purify cDNA Construct

Preparation

- 1 [Optional] Label a new 200 μ l PCR tube 0.1X Pellet Paint.
- 2 Determine the volume of 1X TBE Buffer needed for gel electrophoresis. Dilute the 5X Novex TBE Buffer to 1X.
- 3 Place 6% Novex TBE gel into the gel electrophoresis unit per manufacturer instructions.

Procedure

- 1 Combine the following volumes in the 0.1X Pellet Paint tube. Multiply each volume by the number of libraries being prepared. Make 10% extra reagent if you are preparing multiple libraries.
 - 1X Pellet Paint NF Co-Precipitant (0.2 μ l)
 - Ultrapure water (1.8 μ l)
- 2 Pipette to mix, and then centrifuge briefly.
- 3 Combine 2 μ l CRL and 2 μ l DNA loading dye in a new 1.5 ml microcentrifuge tube.
- 4 Pipette to mix.
- 5 Combine 1 μ l HRL and 1 μ l DNA loading dye in a new 1.5 ml microcentrifuge tube.
- 6 Pipette to mix.
- 7 Combine all amplified cDNA construct (typically 48–50 μ l) and 10 μ l DNA Loading Dye in a new 1.5 ml microcentrifuge tube.
- 8 Pipette to mix.
- 9 Load 2 gel lanes with 2 μ l CRL/loading dye mixture.
- 10 Load 1 gel lane with 2 μ l HRL/loading dye mixture.
- 11 Load 2 gel lanes with 25 μ l each of amplified cDNA construct/loading dye mixture.
- 12 Run the gel for 60 minutes at 145 V or until the blue front dye leaves the gel.
- 13 Remove the gel from the unit.
- 14 Open the cassette according to manufacturer instructions and stain the gel with ethidium bromide in a clean container for 2–3 minutes.
- 15 Place the gel breaker tube into a 2 ml microcentrifuge tube.
- 16 View the gel on a Dark Reader transilluminator or a UV transilluminator.
- 17 Using a razor blade, cut out the bands from the 2 lanes that correspond to the adapter-ligated constructs derived from the 22 nt and 30 nt small RNA fragments.
- 18 Place the band into the 0.5 ml gel breaker tube.
- 19 Centrifuge the nested tubes at $20,000 \times g$ for 2 minutes to move the gel through the holes into the 2 ml tube.
- 20 If you are concentrating the final library, skip the remaining steps and proceed to *Add 300 μ l ultrapure water to the gel debris in the 2 ml tube.* on page 8.
- 21 Add 200 μ l ultrapure water to the gel debris.

- 22 Rotate for at least 2 hours to elute the DNA.
- 23 Transfer the eluate and gel debris to the top of a 5 μm filter.
- 24 Centrifuge at 10 seconds at $600 \times g$.
- 25 Add 300 μl ultrapure water to the gel debris in the 2 ml tube.
- 26 Rotate for at least 2 hours to elute the DNA.
- 27 Transfer the eluate and gel debris to the top of a 5 μm filter.
- 28 Centrifuge at $600 \times g$ for 10 seconds, and then discard the filter.
- 29 Add the following volumes to the eluate:
 - ▶ Glycogen (2 μl)
 - ▶ 3M NaOAc (30 μl)
 - ▶ [Optional] 0.1X Pellet Paint (2 μl)
 - ▶ 100% ethanol (975 μl)
- 30 Centrifuge at $20,000 \times g$ at 20 minutes at 4°C .
- 31 Remove and discard the supernatant. Leave the pellet intact.
- 32 If the pellet becomes loose, centrifuge at $20,000 \times g$ for 2 minutes.
- 33 Wash the pellet with 500 μl 70% ethanol.
- 34 Centrifuge at $20,000 \times g$ for 2 minutes.
- 35 Remove and discard the supernatant. Leave the pellet intact.
- 36 With the lid open, place the tube in a 37°C heat block until the pellet is dry (~7 minutes).
- 37 Resuspend the pellet in 10 μl 10 mM Tris-HCl, pH 8.5.

Check Libraries

- 1 Load 1 μ l resuspended construct on an Agilent Technologies 2100 Bioanalyzer using a DNA-specific chip, such as the DNA 1000 or High Sensitivity DNA.
- 2 Check the size, purity, and concentration of the library.

Normalize Libraries

Procedure

- 1 Normalize library concentration to 2 nM using Tris-HCl 10 mM, pH 8.5.
- 2 For storage, add Tween 20 to the library for a final concentration of 0.1% Tween 20.

SAFE STOPPING POINT

If you are stopping, cap the tube and store at -25°C to -15°C for up to 7 days.

Acronyms

| Acronym | Definition |
|---------|---------------------------|
| cDNA | Complementary DNA |
| CRL | Custom RNA Ladder |
| HML | Ligation Buffer |
| HRL | High Resolution Ladder |
| PCR | Polymerase Chain Reaction |
| PML | PCR Mix |
| RA3 | RNA 3' Adapter |
| RA5 | RNA 5' Adapter |
| RIN | RNA Integrity Number |
| RP1 | RNA PCR Primer |
| RPI | RNA PCR Primer Index |
| RTP | RNA RT Primer |
| STP | Stop Solution |
| UHR | Universal Human Reference |

Notes

Technical Assistance

For technical assistance, contact Illumina Technical Support.

Table 1 Illumina General Contact Information

| | |
|----------------|--------------------------|
| Website | www.illumina.com |
| Email | techsupport@illumina.com |

Table 2 Illumina Customer Support Telephone Numbers

| Region | Contact Number | Region | Contact Number |
|---------------|----------------|-----------------|-----------------|
| North America | 1.800.809.4566 | Japan | 0800.111.5011 |
| Australia | 1.800.775.688 | Netherlands | 0800.0223859 |
| Austria | 0800.296575 | New Zealand | 0800.451.650 |
| Belgium | 0800.81102 | Norway | 800.16836 |
| China | 400.635.9898 | Singapore | 1.800.579.2745 |
| Denmark | 80882346 | Spain | 900.812168 |
| Finland | 0800.918363 | Sweden | 020790181 |
| France | 0800.911850 | Switzerland | 0800.563118 |
| Germany | 0800.180.8994 | Taiwan | 00806651752 |
| Hong Kong | 800960230 | United Kingdom | 0800.917.0041 |
| Ireland | 1.800.812949 | Other countries | +44.1799.534000 |
| Italy | 800.874909 | | |

Safety data sheets (SDSs)—Available on the Illumina website at support.illumina.com/sds.html.

Product documentation—Available for download in PDF from the Illumina website. Go to support.illumina.com, select a product, then select **Documentation & Literature**.



Illumina
5200 Illumina Way
San Diego, California 92122 U.S.A.
+1.800.809.ILMN (4566)
+1.858.202.4566 (outside North America)
techsupport@illumina.com
www.illumina.com