

# TruSight Cardio Sequencing Kit

## Protocol Guide

For Research Use Only. Not for use in diagnostic procedures.

Tagment Genomic DNA	3
Clean Up Tagmented DNA	4
Amplify Tagmented DNA	5
Clean Up Amplified DNA	6
Hybridize Probes	7
Capture Hybridized Probes	8
Perform Second Hybridization	9
Perform Second Capture	10
Clean Up Captured Library	11
Amplify Enriched Library	12
Clean Up Amplified Enriched Library	13
Check Enriched Libraries	14
Acronyms	15
Technical Assistance	17

ILLUMINA PROPRIETARY

illumina®

Document # 15069812 v01  
January 2016

This document and its contents are proprietary to Illumina, Inc. and its affiliates ("Illumina"), and are intended solely for the contractual use of its customer in connection with the use of the product(s) described herein and for no other purpose. This document and its contents shall not be used or distributed for any other purpose and/or otherwise communicated, disclosed, or reproduced in any way whatsoever without the prior written consent of Illumina. Illumina does not convey any license under its patent, trademark, copyright, or common-law rights nor similar rights of any third parties by this document.

The instructions in this document must be strictly and explicitly followed by qualified and properly trained personnel in order to ensure the proper and safe use of the product(s) described herein. All of the contents of this document must be fully read and understood prior to using such product(s).

FAILURE TO COMPLETELY READ AND EXPLICITLY FOLLOW ALL OF THE INSTRUCTIONS CONTAINED HEREIN MAY RESULT IN DAMAGE TO THE PRODUCT(S), INJURY TO PERSONS, INCLUDING TO USERS OR OTHERS, AND DAMAGE TO OTHER PROPERTY.

ILLUMINA DOES NOT ASSUME ANY LIABILITY ARISING OUT OF THE IMPROPER USE OF THE PRODUCT(S) DESCRIBED HEREIN (INCLUDING PARTS THEREOF OR SOFTWARE).

© 2016 Illumina, Inc. All rights reserved.

**Illumina, 24sure, BaseSpace, BeadArray, BlueFish, BlueFuse, BlueGnome, cBot, CSeqPro, CytoChip, DesignStudio, Epicentre, ForenSeq, Genetic Energy, GenomeStudio, GoldenGate, HiScan, HiSeq, HiSeq X, Infinium, iScan, iSelect, MiSeq, MiSeqDx, MiSeq FGx, NeoPrep, NextBio, Nextera, NextSeq, Powered by Illumina, SureMDA, TruGenome, TruSeq, TruSight, Understand Your Genome, UYG, VeraCode, verifi, VeriSeq**, the pumpkin orange color, and the streaming bases design are trademarks of Illumina, Inc. and/or its affiliate(s) in the U.S. and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.

# Tagment Genomic DNA

## Preparation

- 1 Preheat a microheating system with midi plate insert to 58°C.

## Procedure

- 1 Quantify gDNA using a fluorometric method.
- 2 Dilute gDNA in Tris-HCl 10 mM, pH 8.5 to a final volume of 10 µl at 5 ng/µl.
- 3 Add the following items in the order listed to a new midi plate.
  - ▶ Normalized gDNA (10 µl)
  - ▶ TD (25 µl)
  - ▶ TDE1 (15 µl)
- 4 Shake at 1800 rpm for 1 minute.
- 5 Centrifuge at  $280 \times g$  for 1 minute.
- 6 Place on the 58°C microheating system with the lid closed for 10 minutes.
- 7 Add 15 µl ST.
- 8 Shake at 1800 rpm for 1 minute.
- 9 Centrifuge at  $280 \times g$  for 1 minute.
- 10 Incubate at room temperature for 4 minutes.

# Clean Up Tagmented DNA

## Procedure

- 1 Add 65  $\mu$ l SPB.
- 2 Shake at 1800 rpm for 1 minute.
- 3 Incubate at room temperature for 8 minutes.
- 4 Centrifuge at  $280 \times g$  for 1 minute.
- 5 Place on a magnetic stand and wait until the liquid is clear (~2 minutes).
- 6 Remove and discard all supernatant.
- 7 Wash 2 times with 200  $\mu$ l 80% EtOH.
- 8 Use a 20  $\mu$ l pipette to remove residual EtOH.
- 9 Air-dry on the magnetic stand for 10 minutes.
- 10 Remove from the magnetic stand.
- 11 Add 22.5  $\mu$ l RSB.
- 12 Shake at 1800 rpm for 1 minute.
- 13 Incubate at room temperature for 2 minutes.
- 14 Centrifuge at  $280 \times g$  for 1 minute.
- 15 Place on a magnetic stand and wait until the liquid is clear (~2 minutes).
- 16 Transfer 20  $\mu$ l supernatant to a new Hard-Shell PCR plate.

# Amplify Tagmented DNA

## Preparation

- 1 Save the following NLM AMP program on the thermal cycler:
  - ▶ Choose the preheat lid option and set to 100°C
  - ▶ 72°C for 3 minutes
  - ▶ 98°C for 30 seconds
  - ▶ 10 cycles of:
    - ▶ 98°C for 10 seconds
    - ▶ 60°C for 30 seconds
    - ▶ 72°C for 30 seconds
  - ▶ 72°C for 5 minutes
  - ▶ Hold at 10°C

## Procedure

- 1 Arrange Index 1 (i7) adapters in columns 1–12 of the TruSeq Index Plate Fixture.
- 2 Arrange Index 2 (i5) adapters in rows A–H of the TruSeq Index Plate Fixture.
- 3 Place the plate on the TruSeq Index Plate Fixture.
- 4 Using a multichannel pipette, add 5 µl of each Index 1 (i7) adapter down each column. Replace the cap on each i7 adapter tube with a new orange cap.
- 5 Using a multichannel pipette, add 5 µl of each Index 2 (i5) adapter across each row. Replace the cap on each i5 adapter tube with a new white cap.
- 6 Add 20 µl NLM.
- 7 Shake at 1200 rpm for 1 minute.
- 8 Centrifuge at 280 × g for 1 minute.
- 9 Place on the thermal cycler and run the NLM AMP program.

### SAFE STOPPING POINT

If you are stopping, seal the plate and store at 2°C to 8°C for up to 2 days. Alternatively, leave on the thermal cycler overnight.

# Clean Up Amplified DNA

## Procedure

- 1 Centrifuge at  $280 \times g$  for 1 minute.
- 2 Transfer 50  $\mu\text{l}$  supernatant to a new midi plate.
- 3 Add 90  $\mu\text{l}$  SPB.
- 4 Shake at 1800 rpm for 1 minute.
- 5 Incubate at room temperature for 10 minutes.
- 6 Centrifuge at  $280 \times g$  for 1 minute.
- 7 Place on a magnetic stand and wait until the liquid is clear (~2 minutes).
- 8 Remove and discard all supernatant.
- 9 Wash 2 times with 200  $\mu\text{l}$  80% EtOH.
- 10 Use a 20  $\mu\text{l}$  pipette to remove residual EtOH.
- 11 Air-dry on the magnetic stand for 10 minutes.
- 12 Add 27  $\mu\text{l}$  RSB.
- 13 Shake at 1800 rpm for 1 minute.
- 14 Incubate at room temperature for 2 minutes.
- 15 Centrifuge at  $280 \times g$  for 1 minute.
- 16 Place on a magnetic stand and wait until the liquid is clear (~2 minutes).
- 17 Transfer 25  $\mu\text{l}$  supernatant to a new Hard-Shell PCR plate.
- 18 Quantify the library using a fluorometric method.

## SAFE STOPPING POINT

If you are stopping, seal the plate and store at  $-25^{\circ}\text{C}$  to  $-15^{\circ}\text{C}$  for up to 14 days.

# Hybridize Probes

## Preparation

- 1 Save the NRC HYB program on the thermal cycler:
  - ▶ Choose the preheat lid option and set to 100°C
  - ▶ 95°C for 10 minutes
  - ▶ 18 cycles of 1 minute each, starting at 94°C, then decreasing 2°C per cycle
  - ▶ Hold at 58°C

## Pool Libraries

- 1 Combine 500 ng of each DNA library. Make sure that each library has a unique index.
  - ▶ If the total volume is > 40 µl, concentrate the pooled sample to 40 µl.
  - ▶ If the total volume is < 40 µl, increase the volume to 40 µl with RSB.

## Procedure

- 1 Add the following items in the order listed to a new Hard-Shell PCR plate.
  - ▶ DNA library sample or pool (40 µl)
  - ▶ EHB (50 µl)
  - ▶ TCO (10 µl)
- 2 Shake at 1200 rpm for 1 minute.
- 3 Centrifuge at 280 × g for 1 minute.
- 4 Place on the thermal cycler and run the NRC HYB program.
- 5 Keep at the 58°C holding temperature for at least 90 minutes and up to 24 hours.

# Capture Hybridized Probes

## Preparation

- 1 Preheat a microheating system with midi plate insert to 50°C.

## Procedure

- 1 Centrifuge at  $280 \times g$  for 1 minute.
- 2 Transfer all volumes to a new plate.
- 3 Add 250  $\mu$ l SMB.
- 4 Shake at 1200 rpm for 5 minutes.
- 5 Incubate at room temperature for 25 minutes.
- 6 Centrifuge at  $280 \times g$  for 1 minute.
- 7 Place on a magnetic stand and wait until the liquid is clear (~2 minutes).
- 8 Remove and discard all supernatant.
- 9 Remove from the magnetic stand.
- 10 Wash 2 times with 200  $\mu$ l EWS.
- 11 Mix 28.5  $\mu$ l EE1 and 1.5  $\mu$ l HP3, and then vortex.
- 12 Add 23  $\mu$ l elution premix.
- 13 Shake at 1800 rpm for 2 minutes.
- 14 Incubate at room temperature for 2 minutes.
- 15 Centrifuge at  $280 \times g$  for 1 minute.
- 16 Place on a magnetic stand and wait until the liquid is clear (~2 minutes).
- 17 Transfer 21  $\mu$ l supernatant to a new Hard-Shell PCR plate.
- 18 Add 4  $\mu$ l ET2.
- 19 Shake at 1200 rpm for 1 minute.
- 20 Centrifuge at  $280 \times g$  for 1 minute.

## SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to 7 days.



## Perform Second Hybridization

### Procedure

- 1 Add the following reagents in the order listed.
  - ▶ RSB (15  $\mu$ l)
  - ▶ EHB (50  $\mu$ l)
  - ▶ TCO (10  $\mu$ l)
- 2 Shake at 1200 rpm for 1 minute.
- 3 Centrifuge at  $280 \times g$  for 1 minute.
- 4 Place on the thermal cycler and run the NRC HYB program.
- 5 Keep at the 58°C holding temperature for at least 14.5 hours and up to 24 hours.

# Perform Second Capture

## Preparation

- 1 Preheat a microheating system with midi plate insert to 50°C.

## Procedure

- 1 Centrifuge at  $280 \times g$  for 1 minute.
- 2 Transfer supernatant to a new midi plate.
- 3 Add 250  $\mu\text{l}$  SMB.
- 4 Shake at 1200 rpm for 5 minutes.
- 5 Incubate at room temperature for 25 minutes.
- 6 Centrifuge at  $280 \times g$  for 1 minute.
- 7 Place on a magnetic stand and wait until the liquid is clear (~2 minutes).
- 8 Remove and discard all supernatant.
- 9 Remove from the magnetic stand.
- 10 Wash 2 times with 200  $\mu\text{l}$  EWS.
- 11 Mix 28.5  $\mu\text{l}$  EE1 and 1.5  $\mu\text{l}$  HP3, and then vortex.
- 12 Add 23  $\mu\text{l}$  elution premix.
- 13 Shake at 1800 rpm for 2 minutes.
- 14 Incubate at room temperature for 2 minutes.
- 15 Centrifuge at  $280 \times g$  for 1 minute.
- 16 Place on a magnetic stand and wait until the liquid is clear (~2 minutes).
- 17 Transfer 21  $\mu\text{l}$  supernatant to a new midi plate.
- 18 Add 4  $\mu\text{l}$  ET2.
- 19 Shake at 1800 rpm for 1 minute.
- 20 Centrifuge at  $280 \times g$  for 1 minute.

## Clean Up Captured Library

### Procedure

- 1 Add 45  $\mu$ l SPB.
- 2 Shake at 1800 rpm for 1 minute.
- 3 Incubate at room temperature for 10 minutes.
- 4 Centrifuge at  $280 \times g$  for 1 minute.
- 5 Place on a magnetic stand and wait until the liquid is clear (~2 minutes).
- 6 Remove and discard all supernatant.
- 7 Wash 2 times with 200  $\mu$ l 80% EtOH.
- 8 Use a 20  $\mu$ l pipette to remove residual EtOH.
- 9 Air-dry on the magnetic stand for 10 minutes.
- 10 Add 27.5  $\mu$ l RSB.
- 11 Shake at 1800 rpm for 1 minute.
- 12 Incubate at room temperature for 2 minutes.
- 13 Centrifuge at  $280 \times g$  for 1 minute.
- 14 Place on a magnetic stand and wait until the liquid is clear (~2 minutes).
- 15 Transfer 25  $\mu$ l supernatant to a new Hard-Shell PCR plate.

### SAFE STOPPING POINT

If you are stopping, seal the plate and store at  $-25^{\circ}\text{C}$  to  $-15^{\circ}\text{C}$  for up to 7 days.

# Amplify Enriched Library

## Preparation

- 1 Save the following NEM AMP12 program on the thermal cycler:
  - ▶ Choose the preheat lid option and set to 100°C
  - ▶ 98°C for 30 seconds
  - ▶ 12 cycles of:
    - ▶ 98°C for 10 seconds
    - ▶ 60°C for 30 seconds
    - ▶ 72°C for 30 seconds
  - ▶ 72°C for 5 minutes
  - ▶ Hold at 10°C

## Procedure

- 1 Add 5 µl PPC.
- 2 Add 20 µl NEM.
- 3 Shake at 1200 rpm for 1 minute.
- 4 Centrifuge at  $280 \times g$  for 1 minute.
- 5 Place on the thermal cycler and run the NEM AMP12 program.

### **SAFE STOPPING POINT**

If you are stopping, seal the plate and store at 2°C to 8°C for up to 2 days.

## Clean Up Amplified Enriched Library

### Procedure

- 1 Centrifuge at  $280 \times g$  for 1 minute.
- 2 Transfer 50  $\mu\text{l}$  to a new midi plate.
- 3 Add 90  $\mu\text{l}$  SPB.
- 4 Shake at 1800 rpm for 1 minute.
- 5 Incubate at room temperature for 10 minutes.
- 6 Centrifuge at  $280 \times g$  for 1 minute.
- 7 Place on a magnetic stand and wait until the liquid is clear (~2 minutes).
- 8 Remove and discard all supernatant.
- 9 Wash 2 times with 200  $\mu\text{l}$  80% EtOH.
- 10 Use a 20  $\mu\text{l}$  pipette to remove residual EtOH.
- 11 Air-dry on the magnetic stand for 10 minutes.
- 12 Add 32  $\mu\text{l}$  RSB.
- 13 Shake at 1800 rpm for 1 minute.
- 14 Incubate at room temperature for 2 minutes.
- 15 Centrifuge at  $280 \times g$  for 1 minute.
- 16 Place on a magnetic stand and wait until the liquid is clear (~2 minutes).
- 17 Transfer 30  $\mu\text{l}$  supernatant to a new Hard-Shell PCR plate.

### SAFE STOPPING POINT

If you are stopping, seal the plate and store at  $-25^{\circ}\text{C}$  to  $-15^{\circ}\text{C}$  for up to 7 days.

# Check Enriched Libraries

## Quantify Libraries

- 1 Quantify the enriched library using a fluorometric method.

## Assess Quality

- 1 If the library concentration is higher than the supported quantitative range for the High Sensitivity DNA chip, dilute the library 1:10 with RSB.
- 2 Run 1 µl diluted post enriched library using a High Sensitivity DNA chip.
- 3 Denature and dilute pooled libraries to the loading concentration for the instrument you are using. See the denature and dilute libraries guide for your instrument.

## Acronyms

Acronym	Definition
EE1	Enrichment Elution Buffer 1
EHB	Enrichment Hybridization Buffer
ET2	Elute Target Buffer 2
EWS	Enrichment Wash Solution
NEC1	Nextera Enriched Clean Up Plate 1
NEC2	Nextera Enriched Clean Up Plate 2
NEH1	Nextera Enrichment Hyb Plate 1
NEH2	Nextera Enrichment Hyb Plate 2
NEL	Nextera Enrichment Library Plate
NEM	Enrichment Amp Mix
NEW1	Nextera Enrichment Wash Plate 1
NEW2	Nextera Enrichment Wash Plate 2
NIL	Nextera Index Library Plate
NLA	Nextera Library Amplification Plate
NLC	Nextera Library Clean Up Plate
NLM	Library Amp Mix
NLT	Nextera Library Tagment Plate
PPC	PCR Primer Cocktail
RSB	Resuspension Buffer
SMB	Streptavidin Magnetic Beads
SPB	Sample Purification Beads
ST	Stop Tagment Buffer
TCO	TruSight Cardio Oligos
TD	Tagment DNA Buffer
TDE1	Tagment DNA Enzyme TDE

## Notes



## Technical Assistance

For technical assistance, contact Illumina Technical Support.

**Table 1** Illumina General Contact Information

<b>Website</b>	www.illumina.com
<b>Email</b>	techsupport@illumina.com

**Table 2** Illumina Customer Support Telephone Numbers

Region	Contact Number	Region	Contact Number
North America	1.800.809.4566	Japan	0800.111.5011
Australia	1.800.775.688	Netherlands	0800.0223859
Austria	0800.296575	New Zealand	0800.451.650
Belgium	0800.81102	Norway	800.16836
China	400.635.9898	Singapore	1.800.579.2745
Denmark	80882346	Spain	900.812168
Finland	0800.918363	Sweden	020790181
France	0800.911850	Switzerland	0800.563118
Germany	0800.180.8994	Taiwan	00806651752
Hong Kong	800960230	United Kingdom	0800.917.0041
Ireland	1.800.812949	Other countries	+44.1799.534000
Italy	800.874909		

**Safety data sheets (SDSs)**—Available on the Illumina website at [support.illumina.com/sds.html](https://support.illumina.com/sds.html).

**Product documentation**—Available for download in PDF from the Illumina website. Go to [support.illumina.com](https://support.illumina.com), select a product, then select **Documentation & Literature**.



Illumina  
5200 Illumina Way  
San Diego, California 92122 U.S.A.  
+1.800.809.ILMN (4566)  
+1.858.202.4566 (outside North America)  
techsupport@illumina.com  
[www.illumina.com](http://www.illumina.com)