


# TruSight RNA Pan-Cancer Panel

## Protocol Guide

For Research Use Only. Not for use in diagnostic procedures.

Fragment RNA	3
Synthesize First Strand cDNA	4
Synthesize Second Strand cDNA	5
Adenylate 3' Ends	6
Ligate Adapters	7
Perform First PCR Amplification	9
Check Libraries	11
Hybridize Probes	12
Capture Hybridized Probes	13
Perform Second Hybridization	15
Perform Second Capture	16
Clean Up Captured Library	17
Perform Second PCR Amplification	18
Clean Up Amplified Enriched Library	19
Check Enriched Libraries	20
Acronyms	21
Technical Assistance	23



This document and its contents are proprietary to Illumina, Inc. and its affiliates ("Illumina"), and are intended solely for the contractual use of its customer in connection with the use of the product(s) described herein and for no other purpose. This document and its contents shall not be used or distributed for any other purpose and/or otherwise communicated, disclosed, or reproduced in any way whatsoever without the prior written consent of Illumina. Illumina does not convey any license under its patent, trademark, copyright, or common-law rights nor similar rights of any third parties by this document.

The instructions in this document must be strictly and explicitly followed by qualified and properly trained personnel in order to ensure the proper and safe use of the product(s) described herein. All of the contents of this document must be fully read and understood prior to using such product(s).

FAILURE TO COMPLETELY READ AND EXPLICITLY FOLLOW ALL OF THE INSTRUCTIONS CONTAINED HEREIN MAY RESULT IN DAMAGE TO THE PRODUCT(S), INJURY TO PERSONS, INCLUDING TO USERS OR OTHERS, AND DAMAGE TO OTHER PROPERTY.

ILLUMINA DOES NOT ASSUME ANY LIABILITY ARISING OUT OF THE IMPROPER USE OF THE PRODUCT(S) DESCRIBED HEREIN (INCLUDING PARTS THEREOF OR SOFTWARE).

© 2016 Illumina, Inc. All rights reserved.

Illumina, TruSight, the pumpkin orange color, and the streaming bases design are trademarks of Illumina, Inc. and/or its affiliate(s) in the U.S. and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.

Limited Use Label License: This product and its use are the subject of one or more issued and/or pending U.S. and foreign patent applications owned by Max Planck Gesellschaft, exclusively licensed to New England Biolabs, Inc. and sublicensed to Illumina, Inc. The purchase of this product from Illumina, Inc., its affiliates, or its authorized resellers and distributors conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product by the buyer (whether the buyer is an academic or for profit entity). The purchase of this product does not convey a license under any claims in the foregoing patents or patent applications directed to producing the product. The buyer cannot sell or otherwise transfer this product or its components to a third party or otherwise use this product for the following COMMERCIAL PURPOSES: (1) use of the product or its components in manufacturing; or (2) use of the product or its components for therapeutic or prophylactic purposes in humans or animals.

# Fragment RNA

## Preparation

- 1 Save the following Elution 2-Frag-Prime program on the thermal cycler.
  - ▶ Choose the preheat lid option and set to 100°C
  - ▶ 94°C for 8 minutes
  - ▶ Hold at 4°C
  - ▶ Each well contains 17  $\mu$ l
- 2 Set the centrifuge to 15°C to 25°C.

## Procedure

- 1 Dilute the total RNA in nuclease-free ultrapure water to a final volume of 8.5  $\mu$ l in the DFP plate.
- 2 Add 8.5  $\mu$ l EPH.
- 3 Pipette to mix.
- 4 Apply the seal and centrifuge at 280  $\times$  g for 1 minute.



### WARNING

If starting with FFPE RNA, do not perform the following incubation procedure. Proceed immediately to *Synthesize First Strand cDNA* on page 4.

- 5 Place on the thermal cycler and run the Elution 2-Frag-Prime program.
- 6 Centrifuge at 280  $\times$  g for 1 minute.

# Synthesize First Strand cDNA



## WARNING

FSA contains Actinomycin D, a toxin. Personal injury can occur through inhalation, ingestion, skin contact, and eye contact. Dispose of containers and any unused contents in accordance with the governmental safety standards for your region. See the safety data sheet (SDS) for environmental, health, and safety information. For more information, see *Technical Assistance* on page 23.

## Preparation

- 1 Save the following Synthesize 1st Strand program on the thermal cycler:
  - ▶ Choose the preheat lid option and set to 100°C
  - ▶ 25°C for 10 minutes
  - ▶ 42°C for 15 minutes
  - ▶ 70°C for 15 minutes
  - ▶ Hold at 4°C
  - ▶ Each well contains 25  $\mu$ l

## Procedure

- 1 Add 50  $\mu$ l Protoscript II to FSA. Pipette or invert to mix. Then apply the seal and centrifuge briefly.
- 2 Add 8  $\mu$ l Protoscript II and FSA mixture.
- 3 Pipette to mix.
- 4 Place on the thermal cycler and run the Synthesize 1st Strand program.

# Synthesize Second Strand cDNA

## Preparation

- 1 Save the following Synthesize 2nd Strand program on the thermal cycler:
  - ▶ Choose the preheat lid option and set to 30°C
  - ▶ 16°C for 30 minutes
  - ▶ Hold at 4°C
  - ▶ Each well contains 50  $\mu$ l

## Procedure

- 1 Add 5  $\mu$ l RSB.
- 2 Add 20  $\mu$ l SMM.
- 3 Pipette to mix.
- 4 Apply the seal and centrifuge at 280  $\times$  g for 1 minute.
- 5 Place on the preprogrammed thermal cycler and run the Synthesize 2nd Strand program.
- 6 Place on the bench and let stand to bring to room temperature (~5 minutes).
- 7 Add 90  $\mu$ l AMPure XP Beads to the CCP plate.
- 8 Transfer all to the CCP plate.
- 9 Apply the seal and shake at 1800 rpm for 2 minutes.
- 10 Incubate at room temperature for 5 minutes.
- 11 Centrifuge at 280  $\times$  g for 1 minute.
- 12 Place on a magnetic stand and wait until the liquid is clear (~5 minutes).
- 13 Remove and discard 135  $\mu$ l supernatant.
- 14 Wash 2 times with 200  $\mu$ l 80% EtOH.
- 15 Use a 20  $\mu$ l pipette to remove residual EtOH.
- 16 Air-dry on the magnetic stand for 5 minutes.
- 17 Remove from the magnetic stand.
- 18 Add 20  $\mu$ l RSB.
- 19 Apply the seal and shake at 1800 rpm for 2 minutes.
- 20 Incubate at room temperature for 2 minutes.
- 21 Centrifuge at 280  $\times$  g for 1 minute.
- 22 Place on a magnetic stand and wait until the liquid is clear (~5 minutes).
- 23 Transfer 17.5  $\mu$ l supernatant to the ALP plate.

### SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to 7 days.

# Adenylate 3' Ends

## Preparation

- 1 Save the following ATAIL70 program on the thermal cycler:
  - ▶ Choose the preheat lid option and set to 100°C
  - ▶ 37°C for 30 minutes
  - ▶ 70°C for 5 minutes
  - ▶ Hold at 4°C
  - ▶ Each well contains 30  $\mu$ l

## Procedure

- 1 Add 12.5  $\mu$ l ATL.
- 2 Pipette to mix.
- 3 Apply the seal and centrifuge at 280  $\times$  g for 1 minute.
- 4 Place on the thermal cycler and start the program ATAIL70.
- 5 Place on ice for 1 minute or until cooled to 2°C to 8°C.

# Ligate Adapters

## Preparation

- 1 Save the following LIG30 program on the thermal cycler:
  - ▶ Choose the preheat lid option and set to 100°C
  - ▶ 30°C for 10 minutes
  - ▶ Hold at 4°C
  - ▶ Each well contains 37.5  $\mu$ l

## Procedure

- 1 Add the following reagents in the order listed.
  - ▶ RSB (2.5  $\mu$ l)
  - ▶ LIG (2.5  $\mu$ l)
  - ▶ RNA adapters (2.5  $\mu$ l)
- 2 Pipette to mix.
- 3 Apply the seal and centrifuge at  $280 \times g$  for 1 minute.
- 4 Place on the thermal cycler and start the program LIG30.
- 5 Add 5  $\mu$ l STL.
- 6 Pipette to mix.
- 7 Apply the seal and centrifuge at  $280 \times g$  for 1 minute.
- 8 Add 42  $\mu$ l AMPure XP Beads to each well of CAP.
- 9 Transfer entire volume (42  $\mu$ l) from ALP plate to CAP.
- 10 Apply the seal and shake at 1800 rpm for 2 minutes.
- 11 Incubate at room temperature for 5 minutes.
- 12 Centrifuge at  $280 \times g$  for 1 minute.
- 13 Place on a magnetic stand and wait until the liquid is clear (2–5 minutes).
- 14 Remove and discard all supernatant.
- 15 Wash 2 times with 200  $\mu$ l 80% EtOH.
- 16 Use a 20  $\mu$ l pipette to remove residual EtOH.
- 17 Air-dry on the magnetic stand for 5 minutes.
- 18 Remove from the magnetic stand.
- 19 Add 22.5  $\mu$ l RSB.
- 20 Apply the seal and shake at 1800 rpm for 2 minutes.
- 21 Incubate at room temperature for 2 minutes.
- 22 Centrifuge at  $280 \times g$  for 1 minute.
- 23 Place on a magnetic stand and wait until the liquid is clear (2–5 minutes).
- 24 Transfer 20  $\mu$ l supernatant to the PCR plate.

## **SAFE STOPPING POINT**

If you are stopping, seal the plate and store at  $-25^{\circ}\text{C}$  to  $-15^{\circ}\text{C}$  for up to 7 days.



## Perform First PCR Amplification

### Preparation

- 1 Save the following PMM AMP program on the thermal cycler:
  - ▶ Choose the preheat lid option and set to 100°C
  - ▶ 98°C for 30 seconds
  - ▶ 15 cycles of:
    - ▶ 98°C for 10 seconds
    - ▶ 60°C for 30 seconds
    - ▶ 72°C for 30 seconds
  - ▶ 72°C for 5 minutes
  - ▶ Hold at 4°C
  - ▶ Each well contains 50  $\mu$ l

### Procedure

- 1 Place the PCR plate on ice and add 5  $\mu$ l PPC.
- 2 Add 25  $\mu$ l PMM.
- 3 Pipette to mix.
- 4 Place on the thermal cycler and run the PCR program.
- 5 Add 50  $\mu$ l AMPure XP Beads to the PPP plate for each well corresponding to a sample in the PCR plate.
- 6 Apply the seal and centrifuge PCR plate at 280  $\times$  g for 1 minute.
- 7 Transfer the entire volume (50  $\mu$ l) to the PPP plate.
- 8 Apply the seal and shake at 1800 rpm for 2 minutes.
- 9 Incubate at room temperature for 5 minutes.
- 10 Centrifuge at 280  $\times$  g for 1 minute.
- 11 Place on a magnetic stand and wait until the liquid is clear (2–5 minutes).
- 12 Remove and discard all supernatant.
- 13 Wash 2 times with 200  $\mu$ l 80% EtOH.
- 14 Use a 20  $\mu$ l pipette to remove residual EtOH.
- 15 Air-dry on the magnetic stand for 5 minutes.
- 16 Remove from the magnetic stand.
- 17 Add 12.5  $\mu$ l RSB.
- 18 Apply the seal and shake at 1800 rpm for 2 minutes.
- 19 Incubate at room temperature for 2 minutes.
- 20 Centrifuge at 280  $\times$  g for 1 minute.
- 21 Place on a magnetic stand and wait until the liquid is clear (2–5 minutes).
- 22 Transfer 12  $\mu$ l supernatant to the TSP1 plate.

## **SAFE STOPPING POINT**

If you are stopping, seal the plate and store at  $-25^{\circ}\text{C}$  to  $-15^{\circ}\text{C}$  for up to 7 days.

## Check Libraries

### Quantify Library

- 1 If using a Standard Sensitivity NGS Fragment Analysis Kit on an Advanced Analytical Fragment Analyzer, run 2  $\mu$ l undiluted DNA library.
- 2 If using a DNA 1000 chip on an Agilent Technologies 2100 Bioanalyzer, run 1  $\mu$ l undiluted DNA library.
- 3 Check the size and purity of the sample. Expect the final product to be a band at ~250–300 bp.
- 4 Calculate the concentration of the library using a region selection of 160–700 bp.

# Hybridize Probes

## Preparation

- 1 Save the RNA HYB program on the thermal cycler:
  - ▶ Choose the preheat lid option and set to 100°C
  - ▶ 95°C for 10 minutes
  - ▶ 18 cycles of 1 minute each, starting at 94°C, then decreasing 2°C per cycle
  - ▶ 58°C for 90 minutes
  - ▶ Hold at 58°C

## Procedure

- 1 Dilute 200 ng of each library in 10  $\mu$ l RSB.
- 2 Add the following items in the order listed to the RAH1 plate for a final volume of 25  $\mu$ l.
  - ▶ 200 ng library (in 10  $\mu$ l RSB)
  - ▶ CT3 (12.5  $\mu$ l)
  - ▶ RPO (2.5  $\mu$ l)
- 3 Apply the seal and shake at 1200 rpm for 1 minute.
- 4 Centrifuge at 280  $\times$  g for 1 minute.
- 5 Place on the thermal cycler and run the RNA HYB program. Each well contains 25  $\mu$ l.

## Capture Hybridized Probes

### Preparation

- 1 Save the following RNA BIND program on the thermal cycler:
  - ▶ Choose the preheat lid option and set to 100°C
  - ▶ 50°C for 20 minutes
  - ▶ Hold at 50°C

### Procedure

- 1 Centrifuge RAH1 at  $280 \times g$  for 1 minute.
- 2 Add 62.5  $\mu$ l SMB.
- 3 Apply the seal and shake at 1200 rpm for 5 minutes.
- 4 Incubate at room temperature for 25 minutes.
- 5 Centrifuge at  $280 \times g$  for 1 minute.
- 6 Place on a magnetic stand and wait until the liquid is clear (2–5 minutes).
- 7 Remove and discard all supernatant.
- 8 Remove from the magnetic stand.
- 9 Add 50  $\mu$ l EEW.
- 10 Apply the seal and centrifuge at  $280 \times g$  for 10 seconds.
- 11 Pipette to mix.
- 12 Apply the seal and shake at 1800 rpm for 4 minutes.
- 13 Place on the thermal cycler and start the program RNA BIND. Each well contains 52.5  $\mu$ l.
- 14 After 20 minutes, immediately place on a magnetic stand and wait until the liquid is clear (~2 minutes).
- 15 Remove and discard all supernatant.
- 16 Remove from the magnetic stand.
- 17 Repeat steps 9–16 for a total of 2 washes.
- 18 Mix 9.5  $\mu$ l EE1 and 0.5  $\mu$ l HP3, and then vortex.
- 19 Add 10  $\mu$ l elution premix.
- 20 Apply the seal and centrifuge at  $280 \times g$  for 10 seconds.
- 21 Shake at 1800 rpm for 2 minutes.
- 22 Incubate at room temperature for 2 minutes.
- 23 Centrifuge at  $280 \times g$  for 1 minute.
- 24 Place on a magnetic stand and wait until the liquid is clear (~2 minutes).
- 25 Transfer 9  $\mu$ l supernatant to the RAH2 plate.
- 26 Add 1.7  $\mu$ l ET2.

27 Apply the seal and shake at 1200 rpm for 1 minute.

28 Centrifuge at  $280 \times g$  for 1 minute.

### **SAFE STOPPING POINT**

If you are stopping, seal the plate and store at  $-25^{\circ}\text{C}$  to  $-15^{\circ}\text{C}$  for up to 7 days.

## Perform Second Hybridization

### Procedure

- 1 Add the following reagents in the order listed.
  - ▶ CT3 (12.5  $\mu$ l)
  - ▶ RPO (2.5  $\mu$ l)
- 2 Apply the seal and shake at 1200 rpm for 1 minute.
- 3 Centrifuge at  $280 \times g$  for 1 minute.
- 4 Place on the thermal cycler and run the RNA HYB program. Each well contains 25.7  $\mu$ l.

# Perform Second Capture

## Preparation

### Procedure

- 1 Centrifuge RAH2 at  $280 \times g$  for 1 minute.
- 2 Add 62.5  $\mu\text{l}$  SMB.
- 3 Apply the seal and shake at 1200 rpm for 5 minutes.
- 4 Incubate at room temperature for 25 minutes.
- 5 Centrifuge at  $280 \times g$  for 1 minute.
- 6 Place on a magnetic stand and wait until the liquid is clear (2–5 minutes).
- 7 Remove and discard all supernatant.
- 8 Remove from the magnetic stand.
- 9 Add 50  $\mu\text{l}$  EEW.
- 10 Apply the seal and centrifuge at  $280 \times g$  for 10 seconds.
- 11 Pipette to mix.
- 12 Apply the seal and shake at 1800 rpm for 4 minutes.
- 13 Place on the thermal cycler and start the program RNA BIND. Each well contains 53.2  $\mu\text{l}$ .
- 14 After 20 minutes, immediately place on a magnetic stand and wait until the liquid is clear (~2 minutes).
- 15 Remove and discard all supernatant.
- 16 Remove from the magnetic stand.
- 17 Repeat steps 9–16 for a total of 2 washes.
- 18 Mix 9.5  $\mu\text{l}$  EE1 and 0.5  $\mu\text{l}$  HP3, and then vortex.
- 19 Add 10  $\mu\text{l}$  elution premix.
- 20 Apply the seal and centrifuge at  $280 \times g$  for 10 seconds.
- 21 Shake at 1800 rpm for 2 minutes.
- 22 Incubate at room temperature for 2 minutes.
- 23 Centrifuge at  $280 \times g$  for 1 minute.
- 24 Place on a magnetic stand and wait until the liquid is clear (~2 minutes).
- 25 Transfer 9  $\mu\text{l}$  supernatant to the RAW1 plate.
- 26 Add 1.7  $\mu\text{l}$  ET2.
- 27 Apply the seal and shake at 1200 rpm for 1 minute.
- 28 Centrifuge at  $280 \times g$  for 1 minute.



## Clean Up Captured Library

### Preparation

### Procedure

- 1 Add 20  $\mu$ l AMPure XP Beads.
- 2 Apply the seal and shake at 1800 rpm for 2 minutes.
- 3 Incubate at room temperature for 5 minutes.
- 4 Centrifuge at  $280 \times g$  for 1 minute.
- 5 Place on a magnetic stand and wait until the liquid is clear (~2 minutes).
- 6 Remove and discard 27.5  $\mu$ l supernatant.
- 7 Wash 2 times with 200  $\mu$ l 80% EtOH.
- 8 Use a 20  $\mu$ l pipette to remove residual EtOH.
- 9 Air-dry on the magnetic stand for 5 minutes.
- 10 Remove from the magnetic stand.
- 11 Add 27.5  $\mu$ l RSB.
- 12 Apply the seal and shake at 1800 rpm for 2 minutes.
- 13 Incubate at room temperature for 2 minutes.
- 14 Centrifuge at  $280 \times g$  for 1 minute.
- 15 Place on a magnetic stand and wait until the liquid is clear (~2 minutes).
- 16 Transfer 25  $\mu$ l supernatant to the PCR2 plate.

### SAFE STOPPING POINT

If you are stopping, seal the plate and store at  $-25^{\circ}\text{C}$  to  $-15^{\circ}\text{C}$  for up to 7 days.

# Perform Second PCR Amplification

## Preparation

- 1 Save the following EPM AMP program on the thermal cycler:
  - ▶ Choose the preheat lid option and set to 100°C
  - ▶ 98°C for 30 seconds
  - ▶ 14 cycles of:
    - ▶ 98°C for 10 seconds
    - ▶ 60°C for 30 seconds
    - ▶ 72°C for 30 seconds
  - ▶ 72°C for 5 minutes
  - ▶ Hold at 10°C
  - ▶ Each well contains 50  $\mu$ l

## Procedure

- 1 Add 5  $\mu$ l PPC.
- 2 Add 20  $\mu$ l EPM.
- 3 Pipette to mix.
- 4 Apply the seal and centrifuge at 280  $\times$  g for 1 minute.
- 5 Place on the thermal cycler and run the EPM AMP program.

### **SAFE STOPPING POINT**

If you are stopping, seal the plate and store at 2°C to 8°C for up to 2 days.

## Clean Up Amplified Enriched Library

### Preparation

### Procedure

- 1 Centrifuge the PCR2 plate at  $280 \times g$  for 1 minute.
- 2 Add 90  $\mu\text{l}$  AMPure XP Beads to the RAC2 plate.
- 3 Transfer 50  $\mu\text{l}$  from the PCR2 plate to the RAC2 plate.
- 4 Apply the seal and shake RAC2 at 1800 rpm for 2 minutes.
- 5 Incubate at room temperature for 5 minutes.
- 6 Centrifuge at  $280 \times g$  for 1 minute.
- 7 Place on a magnetic stand and wait until the liquid is clear (2–5 minutes).
- 8 Remove and discard 140  $\mu\text{l}$  supernatant.
- 9 Wash 2 times with 200  $\mu\text{l}$  80% EtOH.
- 10 Use a 20  $\mu\text{l}$  pipette to remove residual EtOH.
- 11 Air-dry on the magnetic stand for 5 minutes.
- 12 Remove from the magnetic stand.
- 13 Add 32  $\mu\text{l}$  RSB.
- 14 Apply the seal and shake at 1800 rpm for 1 minute.
- 15 Incubate at room temperature for 2 minutes.
- 16 Centrifuge at  $280 \times g$  for 1 minute.
- 17 Place on a magnetic stand and wait until the liquid is clear (2–5 minutes).
- 18 Transfer 30  $\mu\text{l}$  supernatant to the RAL plate.

### **SAFE STOPPING POINT**

If you are stopping, seal the plate and store at  $-25^{\circ}\text{C}$  to  $-15^{\circ}\text{C}$  for up to 7 days.

## Check Enriched Libraries

### Quantify Libraries

- 1 Quantify the libraries using qPCR according to the Illumina *Sequencing Library qPCR Quantification Guide* (document # 11322363).

### Check Library Quality

- 1 If using a Standard Sensitivity NGS Fragment Analysis Kit on an Advanced Analytical Fragment Analyzer, run 2  $\mu$ l of the postenriched library.
- 2 If using a DNA 1000 Chip, run 1  $\mu$ l of the postenriched library.
- 3 Check the size and purity of the sample. Expect the final product to be a band at ~250–300 bp.
- 4 Check the size of the library for a distribution of DNA fragments with a size range from ~200 bp–1 kb.
- 5 Denature and dilute pooled libraries to the loading concentration for the instrument you are using. For loading recommendations, see the TruSight RNA Pan-Cancer Panel support page.

## Acronyms

Acronym	Definition
ALP	Adapter Ligation Plate
ATL	A-Tailing Mix
CAP	Clean Up ALP Plate
CCP	cDNA Clean Up Plate
CPP	Clean Up PCR Plate
CT3	Capture Target Buffer 3
DFP	Depleted RNA Fragmentation Plate
EE1	Enrichment Elution Buffer 1
EEW	Enhanced Enrichment Wash Buffer
EPH	Elute, Prime, Fragment High Mix
EPM	Enhanced PCR Mix
ET2	Elute Target Buffer 2
FSA	First Strand Synthesis Act D Mix
HP3	2N NaOH
LIG	Ligation Mix
PCR	Polymerase Chain Reaction Plate
PMM	PCR Master Mix
PPC	PCR Primer Cocktail
RAA	RNA Access Amplification Plate
RAC1	RNA Access Clean Up Plate 1
RAC2	RNA Access Clean Up Plate 2
RAH1	RNA Access Hyb Plate 1
RAH2	RNA Access Hyb Plate 2
RAL	RNA Access Library Plate
RAW1	RNA Access Wash Plate 1
RPO	RNA PanCancer Oligos
RSB	Resuspension Buffer
SMB	Streptavidin Magnetic Beads
SMM	Second Strand Marking Master Mix
STL	Stop Ligation Buffer
TSP	Target Sample Plate

## Notes

## Technical Assistance

For technical assistance, contact Illumina Technical Support.

**Table 1** Illumina General Contact Information

<b>Website</b>	www.illumina.com
<b>Email</b>	techsupport@illumina.com

**Table 2** Illumina Customer Support Telephone Numbers

Region	Contact Number	Region	Contact Number
North America	1.800.809.4566	Japan	0800.111.5011
Australia	1.800.775.688	Netherlands	0800.0223859
Austria	0800.296575	New Zealand	0800.451.650
Belgium	0800.81102	Norway	800.16836
China	400.635.9898	Singapore	1.800.579.2745
Denmark	80882346	Spain	900.812168
Finland	0800.918363	Sweden	020790181
France	0800.911850	Switzerland	0800.563118
Germany	0800.180.8994	Taiwan	00806651752
Hong Kong	800960230	United Kingdom	0800.917.0041
Ireland	1.800.812949	Other countries	+44.1799.534000
Italy	800.874909		

**Safety data sheets (SDSs)**—Available on the Illumina website at [support.illumina.com/sds.html](http://support.illumina.com/sds.html).

**Product documentation**—Available for download in PDF from the Illumina website. Go to [support.illumina.com](http://support.illumina.com), select a product, then select **Documentation & Literature**.



Illumina

5200 Illumina Way

San Diego, California 92122 U.S.A.

+1.800.809.ILMN (4566)

+1.858.202.4566 (outside North America)

[techsupport@illumina.com](mailto:techsupport@illumina.com)

[www.illumina.com](http://www.illumina.com)