

Denature RNA

- 1 Dilute 10–100 ng total RNA in nuclease-free ultrapure water to 8.5 μ l.
- 2 Add 8.5 μ l EPH3.
- 3 Pipette 10 times.
- 4 Centrifuge at 280 \times g for 3 seconds.
- 5 Place on the thermal cycler and run the DEN_RNA program.
- 6 Centrifuge at 280 \times g for 10 seconds.

Synthesize First Strand cDNA

- 1 Combine the following volumes to prepare First Strand Synthesis Master Mix.
 - ▶ FSA (9 μ l)
 - ▶ RVT (1 μ l)
- 2 Thoroughly pipette First Strand Synthesis Master Mix.
- 3 Add 8 μ l First Strand Synthesis Master Mix.
- 4 Pipette 10 times.
- 5 Centrifuge at 280 \times g for 10 seconds.
- 6 Place on the thermal cycler and run the FSS program.

Synthesize Second Strand cDNA

- 1 Centrifuge at 280 \times g for 10 seconds.
- 2 Invert SMM to mix, and then centrifuge briefly.
- 3 Add 25 μ l SMM.
- 4 Pipette 10 times.
- 5 Centrifuge at 280 \times g for 10 seconds.
- 6 Place on the thermal cycler and run the SSS program.
- 7 Centrifuge at 280 \times g for 10 seconds.
- 8 Add 90 μ l AMPure XP.
- 9 Shake at 2200 rpm for 1 minute.
- 10 Incubate at room temperature for 5 minutes.
- 11 Centrifuge at 280 \times g for 10 seconds, and then unseal.
- 12 Place on the magnetic stand until liquid is clear.
- 13 Remove and discard supernatant.
- 14 Wash beads as follows.
 - a Add 175 μ l fresh 80% EtOH.
 - b Wait 30 seconds.
 - c Remove and discard supernatant.
- 15 Repeat wash a **second** time.
- 16 Remove residual EtOH.
- 17 Air-dry for 2 minutes.
- 18 Remove from the magnetic stand.
- 19 Add 19.5 μ l RSB.
- 20 Shake at 2700 rpm for 1 minute.
- 21 Incubate at room temperature for 2 minutes.
- 22 Centrifuge at 280 \times g for 10 seconds, and then unseal.
- 23 Place on the magnetic stand until liquid is clear.
- 24 Transfer 17.5 μ l supernatant.

SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to 7 days.

Tagment cDNA

- 1 Centrifuge at 280 × g for 10 seconds.
- 2 Combine the following volumes to prepare Tagmentation Master Mix.
 - ▶ TB1 (11.5 µl)
 - ▶ EBLTL (11.5 µl)
 - ▶ Nuclease-free ultrapure water (14.5 µl)
- 3 Thoroughly vortex the Tagmentation Master Mix.
- 4 Add 32.5 µl Tagmentation Master Mix.
- 5 Pipette thoroughly.
- 6 Place on the thermal cycler and run the TAG program.
- 7 Centrifuge at 280 × g for 10 seconds.
- 8 Incubate at room temperature for 2 minutes.
- 9 Add 10 µl ST2.
- 10 Shake at 2200 rpm for 1 minute.
- 11 Incubate at room temperature for 5 minutes.
- 12 Centrifuge at 280 × g for 10 seconds, and then unseal.
- 13 Place on the magnetic stand until liquid is clear.
- 14 Remove and discard supernatant.
- 15 Wash beads as follows.
 - a Remove from the magnetic stand.
 - b Add 100 µl TWB to each well.
 - c Shake at 2000 rpm for 1 minute.
 - d Centrifuge at 280 × g for 3 seconds.
 - e Place on the magnetic stand until liquid is clear.
 - f Remove and discard supernatant.
- 16 Wash beads a **second** time.
- 17 Wash beads a **third** time, *skipping step f*.
- 18 Combine the following volumes to prepare PCR Master Mix.
 - ▶ EPM (23 µl)

- 19 Thoroughly vortex PCR Master Mix.
- 20 Remove and discard TWB supernatant.
- 21 Remove residual TWB.
- 22 Remove from the magnetic stand.
- 23 Add 40 µl PCR Master Mix.
- 24 Pierce the index adapter plate wells.
- 25 Add 10 µl UDP0XXX.
- 26 Shake at 2000 rpm for 1 minute.
- 27 Centrifuge at 280 × g for 3 seconds.
- 28 Place on the thermal cycler and run the TAG_ PCR program.

SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to 7 days.

Clean Up Library

- 1 Centrifuge at 280 × g for 10 seconds.
- 2 Place on the magnetic stand until liquid is clear.
- 3 Transfer 45 µl supernatant.
- 4 Add 81 µl AMPure XP.
- 5 Shake at 2200 rpm for 1 minute.
- 6 Incubate at room temperature for 5 minutes.
- 7 Centrifuge at 280 × g for 10 seconds, and then unseal.
- 8 Place on the magnetic stand until liquid is clear.
- 9 Remove and discard supernatant.
- 10 Wash beads as follows.
 - a Add 175 µl fresh 80% EtOH.
 - b Wait 30 seconds.
 - c Remove and discard all supernatant.
- 11 Wash beads a **second** time.
- 12 Remove residual EtOH.
- 13 Air-dry on the magnetic stand for 2 minutes.
- 14 Remove from the magnetic stand.
- 15 Add 17 µl RSB.
- 16 Shake at 2700 rpm for 1 minute.
- 17 Incubate at room temperature for 2 minutes.
- 18 Centrifuge at 280 × g for 10 seconds, and then unseal.
- 19 Place on the magnetic stand until liquid is clear.
- 20 Transfer 15 µl supernatant.

SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to 30 days.

Normalize Library

- 1 Analyze 1 µl library with the Qubit dsDNA BR Assay Kit.
- 2 **[Optional]** Analyze 1 µl library with the Agilent 2100 Bioanalyzer System and a DNA 1000 Kit.
- 3 **[Respiratory Virus Panel Libraries]**
 - ▶ For one-plex enrichment, transfer 7.5 µl undiluted library to one well.
 - ▶ For three-plex enrichment, dilute three 200 ng libraries to 2.5 µl each.
- 4 **[All Other Libraries]** Dilute libraries in RSB as follows.
 - ▶ For one-plex enrichment, dilute one 200 ng library to 7.5 µl.
 - ▶ For three-plex enrichment, dilute three 200 ng libraries to 2.5 µl each.
- 5 **[Diluted Libraries]** In one well, combine the 200 ng libraries:

| Number of Libraries | Total Mass (ng) | Total Volume (µl) |
|---------------------|-----------------|-------------------|
| 1 | 200 | 7.5 |
| 3 | 600 | 7.5 |

- ▶ If the total volume is > 7.5 µl, concentrate the pooled sample to 7.5 µl.

Hybridize Probes

- 1 Add the following volumes *in the order listed*.
 - ▶ 200 ng library or 600 ng pool (7.5 µl)
 - ▶ NHB2 (12.5 µl)
 - ▶ Enrichment oligos (2.5 µl)
 - ▶ EHB2 (2.5 µl)
- 2 Pipette 10 times to mix.
- 3 Centrifuge at 280 × g for 3 seconds.
- 4 Place on the thermal cycler and run the HYB program.
- 5 Incubate at 58°C for 90 minutes to 24 hours.

Capture Hybridized Probes

- 1 Combine the following volumes to prepare Elution Master Mix.
 - ▶ EE1 (28.5 µl)
 - ▶ HP3 (1.5 µl)
- 2 Thoroughly pipette Elution Master Mix, and then set aside.
- 3 Centrifuge the PCR plate at 280 × g for 10 seconds.
- 4 Add 62.5 µl SMB.
- 5 Pipette until resuspended.
- 6 Place in the 58°C thermal cycler for 15 minutes.
- 7 Immediately do as follows.
 - a Centrifuge at 280 × g for 10 seconds.
 - b Place on the magnetic stand until liquid is clear.
- 8 Remove and discard supernatant.
- 9 Remove from the magnetic stand.
- 10 Add 50 µl preheated EEW.
- 11 Shake at 2400 rpm for 4 minutes.
- 12 Return unused EEW to the microheating system.
- 13 Return the plate to the 58°C thermal cycler for 5 minutes.
- 14 Immediately do as follows.
 - a Centrifuge at 280 × g for 3 seconds.
 - b Place on the magnetic stand until liquid is clear.
- 15 Remove and discard supernatant.
- 16 Remove from the magnetic stand.
- 17 Add 50 µl preheated EEW.
- 18 Shake at 2000 rpm for 1 minute.
- 19 Return unused EEW to the microheating system.

- 20 Return the plate to the 58°C thermal cycler for 5 minutes.
- 21 Immediately do as follows.
 - a Centrifuge at 280 × g for 3 seconds.
 - b Place on the magnetic stand until liquid is clear.
- 22 Remove and discard supernatant.
- 23 Repeat steps 16–22.
- 24 Remove from the magnetic stand.
- 25 Add 50 µl preheated EEW.
- 26 Shake at 2000 rpm for 1 minute.
- 27 Centrifuge at 280 × g for 3 seconds.
- 28 Transfer 50 µl resuspended bead solution.
- 29 Seal and centrifuge at 280 × g for 3 seconds.
- 30 Return to the 58°C thermal cycler for 5 minutes.
- 31 Immediately place on the magnetic stand until liquid is clear.
- 32 Remove and discard supernatant.
- 33 Remove and discard residual EEW.
- 34 Thoroughly pipette Elution Master Mix.
- 35 Remove from the magnetic stand.
- 36 Add 23 µl Elution Master Mix.
- 37 Shake at 2600 rpm for 1 minute.
- 38 Incubate at room temperature for 2 minutes.
- 39 Centrifuge at 280 × g for 10 seconds, and then unseal.
- 40 Place on the magnetic stand until liquid is clear.
- 41 Transfer 21 µl supernatant.
- 42 Add 4 µl ET2.
- 43 Shake at 2000 rpm for 1 minute.

SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to 7 days.

Amplify Enriched Library

- 1 Centrifuge the sealed plate at 280 × g for 10 seconds.
- 2 Add 5 µl PPC.
- 3 Add 20 µl EPM.
- 4 Shake at 2000 rpm for 1 minute.
- 5 Centrifuge at 280 × g for 10 seconds.
- 6 Place on the thermal cycler and run the AMP program.

Clean Up Enriched Library

- 1 Centrifuge at 280 × g for 10 seconds.
- 2 Add 90 µl AMPure XP.
- 3 Shake at 2200 rpm for 1 minute.
- 4 Incubate at room temperature for 5 minutes.
- 5 Centrifuge at 280 × g for 10 seconds, and then unseal.
- 6 Place on the magnetic stand until liquid is clear.
- 7 Remove and discard supernatant.
- 8 Wash beads as follows.
 - a Add 175 µl fresh 80% EtOH.
 - b Wait 30 seconds.
 - c Remove and discard supernatant.
- 9 Wash beads a **second** time.
- 10 Remove residual EtOH.
- 11 Air-dry on the magnetic stand for 2 minutes.
- 12 Remove from the magnetic stand.
- 13 Add 32 µl RSB.
- 14 Shake at 2600 rpm for 1 minute.
- 15 Incubate at room temperature for 2 minutes.
- 16 Centrifuge at 280 × g for 10 seconds, and then unseal.
- 17 Place on the magnetic stand until liquid is clear.
- 18 Transfer 30 µl supernatant.

SAFE STOPPING POINT

If you are stopping, seal the plate and store at -25°C to -15°C for up to 7 days.

Check Enriched Library

- 1 Check the enriched library:
 - ▶ Analyze 1 µl enriched library with the Qubit dsDNA HS Assay kit.
 - ▶ Analyze 1 µl enriched library with the Agilent 2100 Bioanalyzer System and a DNA 1000 Kit.

Dilute Library to the Starting Concentration

- 1 Obtain the molarity value:
 - ▶ **Bioanalyzer quantification only**—Use the molarity value obtained for the library.
 - ▶ **Bioanalyzer and Qubit quantification**—Calculate molarity value using the average size and concentration.
- 2 Using the molarity value, calculate the volumes of RSB and library needed to dilute libraries to the starting concentration for your system.

| Sequencing System | Starting Concentration (nM) | Final Loading Concentration (pM) |
|-----------------------------|-----------------------------|----------------------------------|
| NextSeq 550 and NextSeq 500 | 20 | 0.8 |
| NovaSeq 6000 | 0.6 | 120 |

- 3 Dilute each library to the starting concentration. Combine 10 µl each diluted library in a tube.
- 4 Follow denature and dilute instructions to dilute libraries.