

Infinium HD

FFPE Restore Protocol

FOR RESEARCH USE ONLY

Infinium HD FFPE Workflow	3
Introduction	4
Preparation	6
Steps	7
Technical Assistance	



ILLUMINA PROPRIETARY
Catalog # WG-901-2004
Part # 15014614 Rev C
September 2011

illumina®

This document and its contents are proprietary to Illumina, Inc. and its affiliates ("Illumina"), and are intended solely for the contractual use of its customer in connection with the use of the product(s) described herein and for no other purpose. This document and its contents shall not be used or distributed for any other purpose and/or otherwise communicated, disclosed, or reproduced in any way whatsoever without the prior written consent of Illumina. Illumina does not convey any license under its patent, trademark, copyright, or common-law rights nor similar rights of any third parties by this document.

The instructions in this document must be strictly and explicitly followed by qualified and properly trained personnel in order to ensure the proper and safe use of the product(s) described herein. All of the contents of this document must be fully read and understood prior to using such product(s).

FAILURE TO COMPLETELY READ AND EXPLICITLY FOLLOW ALL OF THE INSTRUCTIONS CONTAINED HEREIN MAY RESULT IN DAMAGE TO THE PRODUCT(S), INJURY TO PERSONS, INCLUDING TO USERS OR OTHERS, AND DAMAGE TO OTHER PROPERTY.

ILLUMINA DOES NOT ASSUME ANY LIABILITY ARISING OUT OF THE IMPROPER USE OF THE PRODUCT(S) DESCRIBED HEREIN (INCLUDING PARTS THEREOF OR SOFTWARE) OR ANY USE OF SUCH PRODUCT(S) OUTSIDE THE SCOPE OF THE EXPRESS WRITTEN LICENSES OR PERMISSIONS GRANTED BY ILLUMINA IN CONNECTION WITH CUSTOMER'S ACQUISITION OF SUCH PRODUCT(S).

FOR RESEARCH USE ONLY


© 2011 Illumina, Inc. All rights reserved.

Illumina, illuminaDx, BeadArray, BeadXpress, cBot, CPro, DASL, DesignStudio, Eco, GAllx, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, Infinium, iSelect, MiSeq, Nextera, Sentrix, Solexa, TruSeq, VeraCode, the pumpkin orange color, and the Genetic Energy streaming bases design are trademarks or registered trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners.

Infinium HD FFPE Workflow

The following figure illustrates the correct workflows for both Infinium HD FFPE assays (genotyping and methylation). This document describes the FFPE Restoration step of the workflow and assumes that all previous steps in the workflow have been completed.

Figure 1 FFPE Workflow, FFPE Restoration Step



Introduction

The Infinium HD FFPE Restore protocol restores degraded FFPE DNA to a state that is amplifiable by the Infinium HD FFPE genotyping or Infinium HD FFPE methylation whole genome amplification protocol.

DNA samples suitable for DNA Restoration are selected using the Illumina Infinium HD FFPE QC assay.



NOTE

Make sure your centrifuge equipment has buckets deep enough to hold the Zymo-Spin MIDI plate assembly. Some types of centrifuge equipment do not provide buckets deep enough or block access to two stacked plates.

Estimated Time

Hands-on time: ~2 hours

Incubation time: 2 hours

Consumables

Item	Quantity	Storage	Supplied By
Amp Mix Restore Reagent (AMR)	1 tube per 24 FFPE DNA samples	-15° to -25°C	Illumina
Convert Master Mix Reagent (CMM)	1 tube per 24 FFPE DNA samples	-15° to -25°C	Illumina
Elution Restore Buffer Reagent (ERB)	1 tube per 24 FFPE DNA samples	-15° to -25°C	Illumina
Primer Pre Restore Reagent (PPR)	1 tube per 24 FFPE DNA samples	-15° to -25°C	Illumina
Pre Restore (PRS) plate, 96-well 0.8 ml microtiter plate (MIDI)	1 plate	-15° to -25°C	User
Restore (RST) plate, 96-well 0.8 ml microtiter plate (MIDI)	1 plate	-15° to -25°C	User
Zymo™ Purification Kit (ZR-96 DNA Clean & Concentrator™-5)	1 plate per 48 FFPE DNA samples		User
MSA5 plate, 96-well 0.8 ml microtiter plate (MIDI)	1 plate		User
0.1N NaOH	4 µl per well containing sample	2° to 8°C	User
Cap mat	3		User
Troughs	Up to 7 (if using a multi-channel pipette)		User
Foil Adhesive Seal			User
DiH ₂ O	4 µl per well containing sample	-15° to -25°C	User
2.5 L Ice bucket	1		User



NOTE

Thaw all reagents completely at room temperature and allow to equilibrate. Once thawed, gently invert each tube several times to thoroughly mix the reagent. Pulse centrifuge each tube to 280 xg to eliminate bubbles and collect reagent at the bottom of the tube.

Preparation

- ▶ Prepare the wash buffer from the Zymo Purification Kit. Dispense 192 ml 100% ethanol to the 48 ml of DNA Wash Buffer concentrate as described in the Zymo instruction manual.
- ▶ Preheat one heat block to 37°C and one heat block to 95°C and allow the temperature to equilibrate.
- ▶ Apply a PRS barcode label to a new MIDI plate.
- ▶ Thaw the PPR and AMR reagent tubes to room temperature. Vortex to mix the contents. Use 1 tube each per 24 FFPE DNA samples.

Steps

- 1 Prepare the PRS plate by transferring 4 to 8 μl of sample genomic DNA to each well to reach a final DNA amount of at least 100 ng.

**NOTE**

If available, up to 250 ng may be used.

**NOTE**

If you are using bisulfate converted samples for the Infinium HD FFPE methylation assay, take the full elution volume (8 μl) from the bisulfate conversion purification step.

- 2 Prepare the 0.1N NaOH solution:
 - a Dispense 4 μl of 0.1N NaOH to each well containing genomic DNA.
 - b Incubate at room temperature for 10 minutes.
 - 3 Dispense 34 μl of PPR to each sample well in the PRS plate.
 - 4 Dispense 38 μl of AMR to each sample well in the PRS plate; seal the plate with a cap mat.
 - 5 Invert the plate ten times and then centrifuge at 280 $\times g$ for 1 minute.
 - 6 Place the sealed PRS plate on the preheated heat block and close the lid. Incubate the plate at 37°C for 1 hour.
 - 7 At the end of the 1-hour incubation, remove and thaw the ERB and CMM reagent tubes to room temperature. Vortex to mix the contents. Use 1 tube each per 24 FFPE DNA samples.
 - 8 Remove the PRS plate from the heat block and centrifuge at 280 $\times g$ for 1 minute.
- A small icon of a pipette tip with a drop of liquid, indicating a note or important instruction.
- NOTE**
- If you are only using one heat block, set the temperature to 95°C.
- 9 Carefully remove the cap mat.
 - 10 From the Zymo Purification Kit, take the Zymo-Spin™ I-96 Plate and place it on the Collection Plate.
 - 11 Add 7 volumes (560 μl) of Zymo DNA Binding Buffer to each volume of DNA sample in the PRS plate.



NOTE

When using multi-channel pipettes, be careful to avoid cross-contamination between wells.

- 12 Pipette mix 5 times and then transfer sample mixtures to the wells of the Zymo-Spin I-96 Plate mounted on a Collection Plate.
- 13 Centrifuge the Zymo-Spin I-96 Plate at 2250 xg for 2 minutes. Discard the flow-through from the Collection Plate.
- 14 Add 600 μ l of Zymo Wash Buffer (with ethanol added) to each well of the Zymo-Spin I-96 Plate mounted on a Collection Plate.



NOTE

The volume of Zymo Wash Buffer provided in this kit is limited. If you are using multi-channel pipettes and troughs keep all remaining Zymo Wash Buffer for use at step 33.

- 15 Centrifuge the Zymo-Spin I-96 Plate at 2250 xg for 2 minutes. Discard the flow-through from the Collection Plate.
- 16 Apply an RST barcode label to a new MIDI plate.
- 17 Transfer the Zymo-Spin I-96 Plate onto the RST plate and dispense 13 μ l of ERB directly to the column matrix in each well.
- 18 Incubate the plate at room temperature for 5 minutes.
- 19 Centrifuge the Zymo-Spin I-96 RST plate assembly at 2250 xg for 1 minute to elute the DNA. There will be 10 μ l of DNA solution generated.



NOTE

If you are processing less than 96 samples the unused wells can be used for following purification steps.

- 20 Seal the RST plate containing 10 μ l of eluted DNA with a foil adhesive seal. Incubate for 2 minutes at 95°C in the heat block.



NOTE

Do not use a heat sealer to seal the RST plate. This may prevent the use of cap mats for plate sealing later in the protocol. We recommend using foil adhesive seal instead of a heat sealer.

- 21 During the 2-minute incubation, fill the 2.5 L ice bucket with ice.

- 22 Immediately after the 2-minute incubation, transfer the RST plate to the ice for 5 minutes. Press the RST plate into the ice to make sure the bottom of all wells contact the ice.

**NOTE**

If you are only using one heat block, set the temperature to 37°C. Make sure the temperature reaches 37°C before you begin the 1-hour incubation in step 26.

- 23 Leaving the RST plate on ice, remove the foil seal and add 10 µl of CMM to each well. Seal the plate with a new cap mat.
- 24 Vortex the RST plate for 1 minute at 1600 rpm.
- 25 Centrifuge the RST plate for 1 minute at 280 xg.
- 26 Incubate the RST plate at 37°C for 1 hour.
- 27 After the 1-hour incubation, centrifuge the RST plate for 1 minute at 280 xg.
- 28 Take a new Zymo-Spin I-96 Plate and place it on the Collection Plate.

**NOTE**

If you are processing 96 samples use a new plate from Zymo Purification Kit. If you're processing 48 samples or less use the unused wells from the previous purification step.

- 29 Remove the cap mat from the RST plate.
- 30 Add 7 volumes (140 µl) of Zymo DNA Binding Buffer to each volume of DNA sample in the RST plate.

**NOTE**

When using multi-channel pipettes, be careful to avoid cross-contamination between wells.

- 31 Pipette mix 5 times and then transfer sample mixtures to the wells of the Zymo-Spin I-96 Plate mounted on a Collection Plate.
- 32 Centrifuge the Zymo-Spin I-96 Plate at 2250 xg for 2 minutes. Discard the flow-through from the Collection Plate.
- 33 Add 600 µl of Zymo Wash Buffer (with ethanol added) to each well of the Zymo-Spin I-96 Plate mounted on a Collection Plate.
- 34 Centrifuge the Zymo-Spin I-96 Plate at 2250 xg for 2 minutes. Discard the flow-through from the Collection Plate.

- 35 Apply a MSA5 barcode label to a new MIDI plate.
- 36 Transfer the Zymo-Spin I-96 Plate onto the MSA5 plate and dispense 10 μ l of DiH_2O directly to the column matrix in each well.
- 37 Incubate the plate at room temperature for 5 minutes.
- 38 Centrifuge the Zymo-Spin I-96 Plate at 2250 $\times g$ for 1 minute to elute the DNA. There will be 8 μ l of purified DNA solution in the MSA5 plate.



NOTE

The MSA5 plate will be used in the Infinium HD FFPE genotyping or Infinium HD FFPE methylation assay as a DNA plate.

- 39 Do one of the following:
 - Proceed to the Infinium HD FFPE genotyping or Infinium HD FFPE methylation assay documentation to begin processing your restored DNA.
 - Seal the MSA5 plate with adhesive film and store at -20°C if you do not plan to proceed to the Infinium HD FFPE genotyping or Infinium HD FFPE methylation assay immediately.

Technical Assistance

For technical assistance, contact Illumina Customer Support.

Table 1 Illumina General Contact Information

Illumina Website	http://www.illumina.com
Email	techsupport@illumina.com

Table 2 Illumina Customer Support Telephone Numbers

Region	Contact Number	Region	Contact Number
North America	1.800.809.4566	Italy	800.874909
Austria	0800.296575	Netherlands	0800.0223859
Belgium	0800.81102	Norway	800.16836
Denmark	80882346	Spain	900.812168
Finland	0800.918363	Sweden	020790181
France	0800.911850	Switzerland	0800.563118
Germany	0800.180.8994	United Kingdom	0800.917.0041
Ireland	1.800.812949	Other countries	+44.1799.534000

MSDSs

Material safety data sheets (MSDSs) are available on the Illumina website at <http://www.illumina.com/msds>.

Product Documentation

If you require additional product documentation, you can obtain PDFs from the Illumina website if PDFs are available. Go to <http://www.illumina.com/support/documentation.ilmn>. When you click on a link, you will be asked to log in to iCom. After you log in, you can view or save the PDF. To register for an iCom account, please visit <https://icom.illumina.com/Account/Register>.

Illumina, Inc.

9885 Towne Centre Drive

San Diego, CA 92121-1975

+1.800.809.ILMN (4566)

+1.858.202.4566 (outside North America) techsupport@illumina.com

www.illumina.com