

Amplify DNA

- 1 Add DNA into either of the following to create a DNA plate:
 - ▶ Midi plate: 20 µl to each DNA well
 - ▶ TCY plate: 10 µl to each DNA well
- 2 Select **MSA1 Tasks | Make MSA1**.
- 3 Enter the **Number of DNA plates**.
- 4 Place the MA1, MA2, and MSM tubes in the robot tube rack.
- 5 Pour 15 ml NaOH into a trough and place on the robot bed.
- 6 Place DNA and MSA1 plates on robot bed.
- 7 Select **Run**.
- 8 Enter the barcode of each DNA plate.
- 9 Place the DNA plates on the robot bed and select **OK**.
- 10 Vortex the MSA1 plate at 1600 rpm for 1 minute.
- 11 Centrifuge at 280 × g at 22°C for 1 minute.
- 12 Remove the cap mat, place the MSA1 plate on the robot bed, and select **OK**.
- 13 When complete, select **OK**.
- 14 Remove and seal the MSA1 plate.
- 15 Vortex the MSA1 plate at 1600 rpm for 1 minute.
- 16 Centrifuge at 280 × g for 1 minute.

Incubate DNA

- 1 Incubate the MSA1 plate for 20–24 hours at 37°C.

Fragment DNA

- 1 Pulse centrifuge the MSA1 plate at 280 × g.
- 2 Select **MSA1 Tasks | Fragment MSA1**.
- 3 Place the MSA1 plate on the robot bed.
- 4 Place FMS tubes in the robot tube rack.
- 5 Select **Run**.
- 6 When complete, select **OK**.
- 7 Remove the plate and seal with a cap mat.
- 8 Vortex at 1600 rpm for 1 minute.
- 9 Pulse centrifuge at 280 × g.
- 10 Place on the 37°C heat block for 1 hour.

SAFE STOPPING POINT

If you are stopping, seal the plate(s), and store at -25°C to -15°C.

Precipitate DNA

- 1 Select **MSA1 Tasks | Precip MSA1**.
- 2 Pulse centrifuge the sealed plate at 280 × g.
- 3 Place the MSA1 plate on the robot bed.
- 4 Place a half reservoir in the frame, and add PM1 as follows:
 - ▶ For 48 samples, add 1 tube PM1
 - ▶ For 96 samples, add 2 tubes PM1
- 5 Select **Run**.
- 6 Remove the MSA1 plate from the robot bed. Do not select **OK**.
- 7 Vortex at 1600 rpm for 1 minute.
- 8 Incubate on the heat block for 5 minutes.
- 9 Pulse centrifuge at 280 × g for 1 minute.
- 10 Set the centrifuge at 4°C.
- 11 Place the MSA1 plate on the robot bed.
- 12 Select **OK**.
- 13 Remove the MSA1 plate from the robot bed and seal.
- 14 Invert 10 times to mix.
- 15 Incubate at 4°C for 30 minutes.
- 16 Place in the centrifuge.
- 17 Centrifuge at 3000 × g for 20 minutes.
- 18 Remove MSA1 plate.
- 19 Make sure that a blue pellet is present.
- 20 Remove and discard the cap mat.
- 21 Quickly invert the plate and drain the supernatant.
- 22 Firmly tap until all wells are free of liquid.
- 23 Place the plate on a tube rack for 1 hour at room temperature.
- 24 Make sure that a blue pellet is still present.

SAFE STOPPING POINT

If you are stopping, seal the plate(s), and store at -25°C to -15°C.

Resuspend DNA

- 1 Select **MSA1 Tasks | Resuspend MSA1**.
- 2 Place the MSA1 plate on the robot bed.
- 3 Place a quarter reservoir in the frame, and add RA1 as follows:
 - ▶ For 48 samples, add 4.5 ml RA1
 - ▶ For 96 samples, add 9 ml RA1
- 4 Select **Run**.
- 5 Remove the MSA1 plate from the robot deck.
- 6 Apply a foil seal to the MSA1 plate.
- 7 Incubate in the Illumina Hybridization Oven for 1 hour.
- 8 Vortex at 1800 rpm for 1 minute.
- 9 Make sure that the pellets are resuspended.
- 10 Pulse centrifuge at 280 × g.

SAFE STOPPING POINT

If you are stopping, store sealed MSA1 plate(s) at 2°C to 8°C for up to 24 hours. If more than 24 hours, store at -25°C to -15°C.

Store sealed RA1 at -25°C to -15°C. If RA1 will be used the next day, seal it, and store it overnight at 4°C.

Hybridize to BeadChip

- 1 Incubate the MSA1 plate on the heat block for 20 minutes.
- 2 Cool at room temperature for 30 minutes.
- 3 Pulse centrifuge at 280 × g.
- 4 Place the gasket into the hybridization chamber.
- 5 Add 400 µl PB2 into each reservoir.
- 6 Place the hybridization chamber insert into the hybridization chamber.
- 7 Immediately cover the chamber with the lid.
- 8 **[Illumina LIMS] Select Select Infinium HD Super | Confirm for Hyb.**
- 9 **[Illumina LIMS] Scan the barcodes.**
- 10 Remove all BeadChips from packaging.
- 11 Place BeadChips into the robot BeadChip alignment fixtures.
- 12 Place the robot BeadChip alignment fixtures onto the robot deck.
- 13 Pulse centrifuge the MSA1 plate at 280 × g.
- 14 Place the MSA1 plate onto the robot deck.
- 15 Select **Run**.
- 16 Place each robot tip alignment guide on top of each robot BeadChip alignment fixture.
- 17 To start the run, select **OK**.
- 18 When complete, select **OK**.
- 19 Remove the robot BeadChip alignment fixtures.
- 20 Make a record of any sections of BeadChip stripes without complete DNA sample coverage.
- 21 Place each BeadChip in a hybridization chamber insert.

- 22 Place the lid on the chamber and secure with the metal clamps.
- 23 **[Illumina LIMS] Select Infinium HD Super | Prepare Hyb Chamber.**
 - a Scan the barcodes.
- 24 Incubate at 48°C for 16–24 hours.

Prepare for Next Day

- 1 Add 330 ml 100% EtOH to the XC4 bottle.
- 2 Resuspend XC4 by adding 100% EtOH and place the bottle on its side on a rocker until BeadChips are ready for coating. Alternatively, leave the bottle upright on the lab bench overnight.
- 3 Soak the tip guide inserts in a 1% aqueous Alconox solution.
- 4 Rinse the tip guides with DiH₂O at least 3 times.
- 5 Dry the tip guide.

Wash BeadChips

This step prepares the BeadChips for the staining process.

- 1 Submerge the wash rack in the PB1 wash.
- 2 Remove the hybridization chamber inserts.
- 3 Remove the BeadChips .
- 4 Remove the cover seals from the BeadChips.
- 5 Place the BeadChips into the submerged wash rack.
- 6 Move the wash rack up and down for 1 minute.
- 7 Move the wash rack to the next PB1Wash.
- 8 Move the wash rack up and down for 1 minute.
- 9 Confirm that you are using the correct Infinium glass back plates and spacers.
- 10 Fill the BeadChip alignment fixture with 150 ml PB1 for up to 8 BeadChips.
- 11 For each BeadChip, place one black frame into the BeadChip alignment fixture.
- 12 Place each BeadChip into a black frame.
- 13 Place a **clear** spacer onto the top of each BeadChip.
- 14 Place the alignment bar onto the alignment fixture.
- 15 Place a clean glass back plate on top of each clear spacer.
- 16 Secure each flow-through chamber assembly with metal clamps.
- 17 Remove the assembled flow-through chamber from the alignment fixture.
- 18 Trim the spacers from each end of the assembly.

- 19 Leave assembled flow-through chambers on the lab bench.
- 20 Wash the hybridization chamber reservoirs with DI H₂O.

Extend and Stain BeadChips

- 1 Select **XStain Tasks | XStain HD BeadChip**.
- 2 If imaging immediately after staining, turn on the scanner.
- 3 Add the following reagents to reservoirs:

Reagent	# BeadChips	Volume
95% formamide/1 mM EDTA	1–8	15 ml
	9–16	17 ml
	17–24	25 ml
RA1	1–8	10 ml
	9–16	20 ml
	17–24	30 ml
XC3	1–8	50 ml
	9–16	100 ml
	17–24	150 ml

- 4 Invert the XC1, XC2, TEM, STM, and ATM tubes to mix. Remove the caps, and place on the robot deck.
- 5 Enter the number of BeadChips.
- 6 Select **Run**.
- 7 **[Non-Illumina LIMS]** Enter the stain temperature listed on the STM tube.
- 8 Place the flow-through chambers into the chamber rack.
- 9 Select **OK**.
- 10 Remove the flow-through chambers from the chamber rack.
- 11 Fill the water circulator.
- 12 Turn on the water circulator and set the temperature to 44°C.

- 13 Set up two top-loading wash dishes labeled PB1 and XC4.
- 14 Add 310 ml PB1 to the PB1 wash dish.
- 15 Submerge the staining rack in the wash dish.
- 16 Leave the staining rack in the wash dish.
- 17 Disassemble each flow-through chamber.
- 18 Place the BeadChips into the submerged staining rack.
- 19 Slowly move the staining rack up and down 10 times.
- 20 Soak for 5 minutes.
- 21 Vigorously shake the XC4 bottle.
- 22 Add 310 ml XC4 to the XC4 wash dish and cover.
- 23 Transfer the staining rack to the XC4 wash dish.
- 24 Slowly lift the staining rack up and down 10 times.
- 25 Soak for 5 minutes.
- 26 Remove the staining rack and place it onto the tube rack.
- 27 Place the tube rack into the vacuum desiccator.
- 28 Dry the BeadChips for 50–55 minutes at 675 mm Hg (0.9 bar).

Acronyms

Acronym	Definition
EDTA	Ethylenediaminetetraacetic acid
EtOH	Ethanol
ATM	Anti-Stain Two-Color Master Mix
FMS	Fragmentation solution
MA1	Multi-Sample Amplification 1 Mix
MA2	Multi-Sample Amplification 2 Mix
MSM	Multi-Sample Amplification Master Mix
PB1	Reagent used to prepare BeadChips for hybridization
PB2	Humidifying buffer used during hybridization
PM1	Precipitation solution
RA1	Resuspension, hybridization, and wash solution
STM	Superior Two-Color Master Mix
TEM	Two-Color Extension Master Mix
XC1	XStain BeadChip solution 1
XC2	XStain BeadChip solution 2
XC3	XStain BeadChip solution 3
XC4	XStain BeadChip solution 4

SAFE STOPPING POINT

Store the BeadChips in the Illumina BeadChip Slide Storage Box at room temperature. Scan within 72 hours.