

# **cBot Software Release Notes**

## **cBot v3.0.46**

### **Recipe Installer v3.0.4**

***For cBot instruments used with HiSeq 2500,  
HiSeq 3000/4000, and HiSeq X systems***

**May 3, 2016**

## Introduction

These release notes describe key changes to software components for cBot 3.0.46 since the release of cBot v3.0.40.

This software release supports the release of the cBot 2 system, which has barcode scanners, and clustering on cBot or cBot 2 using the HiSeq 3000/4000 SR Cluster Kit.

The software package includes:

- cBot v3.0.46
- Recipe Installer 3.0.4

## I. cBot 3.0.46

### NEW FEATURES:

- Added recipe for the HiSeq 3000/4000 single-read (SR) flow cell
- For LIMS integration, the task name is now specified between rehyb and cluster generation.
- Added robustness to software updates by backing up config files.

### DEFECT REPAIRS:

- None

### KNOWN ISSUES:

- The software does not verify whether the flow cell version matches the reagent version.
- When scanning a HiSeq flow cell package to record the flow cell serial number, the software does not make sure that the correct 9-character flow cell serial number is scanned. Make sure that the 9-character serial number is scanned, and not a lot number or other package barcode.
- If a cBot \*.msi file exists in the temp folder before installation (due to a previous installation attempt), the software installer deletes the existing file. If 2 \*.msi files are present, the installer does not proceed as a safeguard preventing installation of the incorrect \*.msi file. Relaunch the installer to complete the installation process.

## II. Recipe Installer 3.0.4

### NEW FEATURES:

- Added the following recipe:
  - HiSeq\_3000\_4000\_SR\_HD\_Exclusion\_Amp\_v1.0.
  
- The following recipes have reinstated since the v2.0.34 release:
  - PE\_HiSeq\_Cluster\_Kit\_v4\_cBot\_Rehyb\_recipe\_v9.0.xml
  - SR\_HiSeq\_Cluster\_Kit\_v4\_cBot\_Rehyb\_recipe\_v9.0.xml