

Master Lot Sheet

Description

Product Name	VeriSeq NIPT SMP Prep Kit(96 SMP)		
Part Number	15066802	Catalog Number	15066802
Lot Number	A117200	Kit Expiration Date	31-JUL-2019

Master Kit Contents

Box 1 of 4

Item Description	Part Number	Lot Number
VeriSeq™ NIPT Extraction Box (96)	15066807	A117200-1
cfDNA Extraction Box (96)	15068028	0020093240

Box 2 of 4

Item Description	Part Number	Lot Number
VeriSeq™ NIPT Library Prep Box(96)	15066810	A117200-2
DX#-VER1, End Repair Mix	20005843	20314200
DX#-VAT1, A-Tailing Mix	20005845	20331089
DX#-VLG1, Ligation Mix	20005847	20331093
DX#-HT1, Hybridization Buffer	15072315	20331342
VeriSeq™ NIPT DNA Adapter Plate	15068029	V164386053

Box 3 of 4

Item Description	Part Number	Lot Number
VeriSeq™ NIPT Accessory Box	15066811	A117200-3
DNA Binding Plate	15072439	1804/002
DX#-RSB, Resuspension Buffer	15067330	20258708
DX#-SPB, Sample Purification Beads	15067744	20331099
DX#-DQR, DNA Quantitation RGT	15067761	20331340
DX#-DQS, DNA Quantitation Standard	15067764	20335694

Box 4 of 4

Item Description	Part Number	Lot Number
VeriSeq™ NIPT Workflow Tubes and Labels	15071543	A117200-4
Empty Pooling Tube	15071858	20330347