

# **Local Run Manager**

# **RNA Amplicon Analysis Module v1.0.2**

## **Release Notes**

**08-JAN-2018**

## Introduction

These Release Notes detail notable items for the Local Run Manager RNA Amplicon Analysis Module v1.0 release.

Please note that the RNA Amplicon Analysis Module v1.0.2 requires Local Run Manager Framework v2.0.0 or higher and is not compatible with Local Run Manager Framework v1.3.1 or lower.

For more information about this analysis module and how to use it, refer to the Local Run Manager RNA Amplicon Analysis Module Workflow Guide, available from the Local Run Manager Support Page, Documentation and Literature, on [illumina.com](http://illumina.com).

[http://support.illumina.com/sequencing/sequencing\\_software/local-run-manager/documentation.html](http://support.illumina.com/sequencing/sequencing_software/local-run-manager/documentation.html)

### NEW FEATURES:

- Initial Release

### DEFECT REPAIRS:

- None

### KNOWN ISSUES:

- In rare instances, the Sample Sheet file generated during analysis may not reflect the same ordering of samples as what was entered on the "Run Setup" page.